SUPPORTING INFECTIOUS DISEASE RESEARCH

# *Campylobacter jejuni* subsp*. jejuni*, Strain MK 104

## Catalog No. NR-4084

(Derived from ATCC<sup>®</sup> 43446<sup>™</sup>)

## For research use only. Not for human use.

Contributor: ATCC<sup>®</sup>

#### Manufacturer:

**BEI Resources** 

#### **Product Description:**

Bacteria Classification: Campylobacteraceae,

Campylobacter

<u>Species</u>: *Campylobacter jejuni* subsp. *jejuni* Strain: MK 104

Serotype: 0:19

- Serotype: 0:19
- <u>Original Source</u>: *Campylobacter jejuni* (*C. jejuni*) subsp. *jejuni*, strain MK 104 was isolated by Dr. M. Karmali from human feces at The Hospital for Sick Children in Toronto, Ontario, Canada.
- <u>Comment</u>: This strain was deposited to the ATCC<sup>®</sup> by Dr. J. L. Penner in 1986.<sup>1</sup>

*C. jejuni* is a Gram-negative slender, curved, motile rod commonly found in animal feces. It is a microaerophilic organism that is very sensitive to environmental stresses.<sup>2</sup> *C. jejuni* is among the most frequently identified bacterial causes of human gastroenteritis in the United States and other industrialized countries.<sup>3</sup> Food poisoning caused by *C. jejuni* can be largely attributed to the consumption of contaminated food animal products, especially poultry. In most cases, the resulting infection can be severely debilitating but is rarely life-threatening. However, in some cases, *C. jejuni* infections have been linked to the subsequent development of two neuropathies, Guillain-Barré syndrome<sup>2,4,5</sup> and Miller-Fisher syndrome<sup>4</sup> and to a reactive arthropathy, Reiter syndrome.<sup>2</sup>

## **Material Provided:**

Each vial contains approximately 0.5 mL of bacterial culture in Brucella broth supplemented with 10% glycerol.

<u>Note</u>: If homogeneity is required for your intended use, please purify prior to initiating work.

## Packaging/Storage:

NR-4084 was packaged aseptically, in screw-capped plastic cryovials. The product is provided frozen and should be stored at -60°C or colder immediately upon arrival. For long-term storage, the vapor phase of a liquid nitrogen freezer is recommended. Freeze-thaw cycles should be avoided.

## **Growth Conditions:**

#### Media:

Brucella broth, Tryptic soy broth or equivalent

Tryptic Soy agar with 5% defibrinated sheep blood or Brucella agar or equivalent

Incubation:

Temperature: 37°C to 42°C

Atmosphere: Microaerophilic (3 to 5% O<sub>2</sub> and 4 to 8% CO<sub>2</sub>) Propagation:

- 1. Keep vial frozen until ready for use, then thaw.
- 2. Transfer the entire thawed aliquot into a single tube of broth.
- 3. Use the suspension to inoculate an agar slant and/or plate.
- 4. Incubate the tube, slant and/or plate at 37°C to 42°C for 24 to 48 hours.
- <u>Note</u>: The thawed vial may be plated directly on TSA with 5% defibrinated sheep blood and grown at 37°C to 42°C in a microaerophilic atmosphere. This may require a longer incubation time than the biphasic culture.

#### Citation:

Acknowledgment for publications should read "The following reagent was obtained through BEI Resources, NIAID, NIH: *Campylobacter jejuni* subsp. *jejuni*, Strain MK 104, NR-4084."

#### **Biosafety Level: 2**

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the following publication: U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. <u>Biosafety in</u> <u>Microbiological and Biomedical Laboratories</u>. 5th ed. Washington, DC: U.S. Government Printing Office, 2009; see <u>www.cdc.gov/biosafety/publications/bmbl5/index.htm</u>.

#### **Disclaimers:**

You are authorized to use this product for research use only. It is not intended for human use.

Use of this product is subject to the terms and conditions of the BEI Resources Material Transfer Agreement (MTA). The MTA is available on our Web site at <u>www.beiresources.org</u>.

While BEI Resources uses reasonable efforts to include accurate and up-to-date information on this product sheet, neither ATCC<sup>®</sup> nor the U.S. Government makes any warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. Neither ATCC<sup>®</sup> nor the U.S. Government warrants that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, use and disposal.  $ATCC^{\textcircled{B}}$  and the U.S. Government are not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to ensure

E-mail: <u>contact@beiresources.org</u> Tel: 800-359-7370 Fax: 703-365-2898 di**c**i**i** resources

SUPPORTING INFECTIOUS DISEASE RESEARCH

authenticity and reliability of materials on deposit, the U.S. Government,  $ATCC^{\textcircled{m}}$ , their suppliers and contributors to BEI Resources are not liable for damages arising from the misidentification or misrepresentation of products.

#### **Use Restrictions:**

This material is distributed for internal research, noncommercial purposes only. This material, its product or its derivatives may not be distributed to third parties. Except as performed under a U.S. Government contract, individuals contemplating commercial use of the material, its products or its derivatives must contact the contributor to determine if a license is required. U.S. Government contractors may need a license before first commercial sale.

#### **References:**

- Mills, S. D., L. A. Kurjanczyk, and J. L. Penner. "Identification of an Antigen Common to Different Species of the Genus *Campylobacter*." <u>J. Clin. Microbiol.</u> 26 (1988): 1411–1413. PubMed: 2457601.
- Altekruse, S. F., et al. "Campylobacter jejuni-An Emerging Foodborne Pathogen." <u>Emerg. Infect. Dis.</u> 5 (1999): 28-35. PubMed: 10081669.
- Gibreel, A. and D. E. Taylor. "Macrolide Resistance in Campylobacter jejuni and Campylobacter coli." <u>J.</u> <u>Antimicrob. Chemother.</u> 58 (2006): 243-255. PubMed: 16735431.
- Woodward, D. L. and F. G. Rodgers. "Identification of *Campylobacter* Heat-Stable and Heat-Labile Antigens by Combining Penner and Lior Serotyping Schemes." <u>J.</u> <u>Clin. Microbiol.</u> 40 (2002): 741–745. PubMed: 11880386.
- Sinha, S., et al. "Detection of Preceding Campylobacter jejuni Infection by Polymerase Chain Reaction in Patients with Guillain-Barre Syndrome." <u>Trans. R. Soc.</u> <u>Trop. Med. Hyg.</u> 98 (2004): 342-346. PubMed: 15099989.

ATCC<sup>®</sup> is a trademark of the American Type Culture Collection.

