

Certificate of Analysis for NR-46220

Staphylococcus aureus, Strain GA-62

Catalog No. NR-46220

Product Description: Staphylococcus aureus (S. aureus), strain GA-62 was isolated in 2005 from the blood of a 72 year-old-male with cellulitis and/or a bloodstream infection in Georgia, USA. S. aureus, strain GA-62 is a methicillin-resistant S. aureus (MRSA) strain.

Lot¹: 63693265 Manufacturing Date: 21AUG2015

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-positive cocci	Gram-positive cocci
Colony morphology ²	Report results	Circular, flat, entire, smooth and
Colony morphology	Report results	gray (Figure 1)
Motility (wet mount)	Report results	Non-motile
Hemolysis ²	Report results	β-hemolytic
Biochemical characterization	Report results	p-nemolytic
Catalase	Positive	Positive
Catalase Coagulase ³		
VITEK [®] 2 Compact (GP card)	Report results Consistent with S. aureus	Positive
VITER 2 Compact (GP card)	Consistent with S. aureus	Consistent with S. aureus
Antibiotic Susceptibility Profile		
VITEK [®] (AST-GP71 card) ⁴		
Beta-lactamase ⁵	Report results	Positive
Cefoxitin screen	Report results	Positive
Benzylpenicillin	Report results	Resistant (≥ 0.5 µg/mL)
Oxacillin	Resistant	Resistant (≥ 4 µg/mL)
Gentamicin	Sensitive	Sensitive (≤ 0.5 µg/mL)
Ciprofloxacin	Report results	Resistant (≥ 8 µg/mL)
Levofloxacin	Resistant	Resistant (≥ 8 µg/mL)
Moxifloxacin	Report results	Resistant (≥ 8 µg/mL)
Clindamycin (inducible resistance)	Report results	Positive ⁶
Erythromycin	Resistant	Resistant (≥ 8 µg/mL)
Quinupristin/dalfopristin	Report results	Sensitive (≤ 0.25 µg/mL)
Linezolid	Sensitive	Sensitive (= 2 µg/mL)
Daptomycin	Sensitive	Sensitive (= 1 µg/mL)
Vancomycin	Sensitive	Sensitive (= 1 µg/mL)
Minocycline	Report results	Sensitive (≤ 0.5 µg/mL)
Tetracycline	Sensitive	Sensitive (≤ 1 μg/mL)
Tigecycline	Report results	Sensitive (≤ 0.12 μg/mL) ⁷
Nitrofurantoin	Report results	Sensitive (≤ 16 µg/mL)
Rifampicin	Sensitive	Sensitive (≤ 0.5 μg/mL)
Trimethoprim/sulfamethoxazole	Resistant	Resistant (≥ 320 µg/mL)
Etest [®] antibiotic test strips ⁸	. toolotant	1.00.0.0 (= 020 pg/112)
Chloramphenicol ⁹	Sensitive	Sensitive (= 4 µg/mL)
Teicoplanin ⁹	Report results	Sensitive (= 4 µg/mL)
·	Roport roduito	σοποιίνο (= 1.5 μg/πιΣ)
Genotypic Analysis		
Sequencing of 16S ribosomal RNA gene	Consistent with S. aureus	Consistent with S. aureus
(~ 1390 base pairs)		
Purity (post-freeze) ¹⁰	Consistent with expected colony	Consistent with expected colony
	morphology	morphology
Viability (post-freeze) ²	Growth	Growth

BEI Resources www.beiresources.org E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898



SUPPORTING INFECTIOUS DISEASE RESEARCH

Certificate of Analysis for NR-46220

¹S. aureus, strain GA-62 was deposited to BEI Resources as part of the NARSA collection. NR-46220 was produced by inoculation of the deposited material into Tryptic Soy broth and grown 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar with 5% defibrinated sheep blood kolles which were grown 1 day at 37°C in an aerobic atmosphere to produce this lot.

²1 day at 37°C in an aerobic atmosphere on Tryptic Soy agar with 5% defibrinated sheep blood

³4 hours at 37°C in rabbit serum with 0.15% EDTA (Coagulase Plasma BBL™ 240827)

⁴Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S22 (2012)

⁵The production of beta-lactamase was detected using a Cefinase™ Paper Disc (BBL™ 231650).

⁶The VITEK® AST-GP71 card tests for both clindamycin resistance and inducible clindamycin resistance (ICR). A positive ICR test is indicative of inducible MLS_b resistance, which confers resistance to macrolides, lincosamides, and type B streptogramin and the isolate should be considered resistant to clindamycin. *S. aureus*, strain GA-62 was found to be sensitive to clindamycin but had a positive ICR test and therefore is considered resistant to clindamycin.


⁷MIC Interpretation Guideline: EUCAST Version 4.0 (2014)

81 day at 37°C in an aerobic atmosphere on Mueller Hinton agar

⁹For both chloramphenicol (bioMérieux Etest[®] 412308) and teicoplanin (bioMérieux Etest[®] 412459), a MIC ≤ 8 μg/mL is sensitive, a MIC = 16 μg/mL is intermediate and a MIC ≥ 32 μg/mL is resistant.

¹⁰Purity of this lot was assessed for 7 days at 37°C in an aerobic atmosphere on Tryptic Soy agar with 5% defibrinated sheep blood.

Figure 1: Colony Morphology



Date: 30 NOV 2015

Signature:

BEI Resources Authentication

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898