

Certificate of Analysis for NR-46253

Staphylococcus aureus, Strain OR-229

Catalog No. NR-46253

Product Description: *Staphylococcus aureus (S. aureus*), strain OR-229 is of unknown origin. *S. aureus*, strain OR-229 is a clinically-associated methicillin-resistant *S. aureus* (MRSA) strain.

Lot¹: 64044908 Manufacturing Date: 02MAR2016

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-positive cocci	Gram-positive cocci
Colony morphology ²	Report results	Circular, convex, entire, smooth and
Colony morphology	1 toport rooms	yellow (Figure 1)
Motility (wet mount)	Report results	Non-motile
Hemolysis ³	Report results	β-hemolytic
Biochemical characterization	Roport rodato	p nemery ac
Catalase	Positive	Positive
Coagulase ⁴	Report results	Positive
VITEK® 2 Compact (GP card)	≥ 90% probability of being S. aureus	S. aureus (99% probability) ⁵
VITER 2 Compact (Gr card)	2 30 % probability of being 6. aureus	G. dureus (3976 probability)
Antibiotic Susceptibility Profile		
VITEK® (AST-GP71 card)6		
Beta-lactamase ⁷	Report results	Positive
Cefoxitin screen	Report results	Positive
Benzylpenicillin	Report results	Resistant (≥ 0.5 µg/mL)
Oxacillin	Resistant	Resistant (≥ 4 µg/mL)
Gentamicin	Sensitive	Sensitive (≤ 0.5 µg/mL)
Ciprofloxacin	Report results	Resistant (≥ 8 µg/mL)
Levofloxacin	Resistant	Resistant (= 4 µg/mL)
Moxifloxacin	Report results	Intermediate (= 1 µg/mL)
Clindamycin (inducible resistance)	Report results	Negative
Erythromycin	Resistant	Resistant (≥ 8 µg/mL)
Clindamycin	Sensitive	Sensitive (≤ 0.25 µg/mL)
Quinupristin/dalfopristin	Report results	Sensitive (≤ 0.25 µg/mL)
Linezolid	Sensitive	Sensitive (= 2 µg/mL)
Daptomycin	Sensitive	Sensitive (= 0.5 µg/mL)
Vancomycin	Sensitive	Sensitive (≤ 1 µg/mL)
Minocycline	Report results	Sensitive (≤ 0.5 µg/mL)
Tetracycline	Sensitive	Sensitive (≤ 1 µg/mL)
Tigecycline	Report results	Sensitive (≤ 0.12 µg/mL) ⁸
Nitrofurantoin	Report results	Sensitive (≤ 16 µg/mL)
Rifampicin	Sensitive	Sensitive (≤ 0.5 µg/mL)
Trimethoprim/sulfamethoxazole	Sensitive	Sensitive (≤ 10 µg/mL)
Etest® antibiotic test strips9	5.101.10	(= 10 µg/111L)
Chloramphenicol ¹⁰	Sensitive	Sensitive (= 3 µg/mL)
Teicoplanin ¹⁰	Report results	Sensitive (= 0.38 µg/mL)
Totopianii	report results	σοποιανο (= σ.σο μg/π.ε)
Genotypic Analysis	> 000/	00.00/
Sequencing of 16S ribosomal RNA gene	≥ 99% sequence identity to	99.9% sequence identity to
(~ 870 base pairs)	S. aureus type strain	S. aureus type strain
((GenBank: L37597)	(GenBank: L37597)
Descrites (reach frages)11	Consistent with expected colony	Consistent with expected colony
Purity (post-freeze) ¹¹	morphology	morphology
Viability (next fraces)?		
Viability (post-freeze) ²	Growth	Growth

BEI Resources www.beiresources.org E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898



Certificate of Analysis for NR-46253

SUPPORTING INFECTIOUS DISEASE RESEARCH

1S. aureus, strain OR-229 was deposited to BEI Resources as part of the NARSA collection. NR-46253 was produced by inoculation of the deposited material into Tryptic Soy broth and grown 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles which were grown 1 day at 37°C in an aerobic atmosphere to produce this lot.

²1 day at 37°C in an aerobic atmosphere on Tryptic Soy agar

³1 day at 37°C in an aerobic atmosphere on Tryptic Soy agar with 5% defibrinated sheep blood

⁴1 day at 37°C in rabbit serum with 0.15% EDTA (Coagulase Plasma BBL™ 240827)

⁵Percent probabilities above 90% indicate a close match to the typical biochemical pattern for the given organism, with a percent probability of 99% being a perfect match between the test reaction pattern and the unique biochemical pattern of the given organism or organism group. For additional information, please refer to O'Hara, C.M. and J. M. Miller. "Evaluation of the VITEK 2 ID-GNB Assay for Identification of Members of the Family *Enterobacteriaceae* and Other Nonenteric Gram-Negative Bacilli and Comparison with the VITEK GNI+ Card." J. Clin. Microbiol. 41 (2003): 2096-2101. PubMed: 12734254.

⁶Minimum Inhibitory Concentration (MIC): MIC Interpretation Guideline: CLSI M100-S22 (2012)

⁷The production of beta-lactamase was detected using a Cefinase™ Paper Disc (BBL™ 231650).

⁸MIC Interpretation Guideline: EUCAST Version 4.0 (2014)

⁹1 day at 37°C in an aerobic atmosphere on Mueller Hinton agar

¹ºFor both chloramphenicol (bioMérieux Etest® 412308) and teicoplanin (bioMérieux Etest® 412459), a MIC ≤ 8 µg/mL is sensitive, a MIC = 16 µg/mL is intermediate and a MIC ≥ 32 µg/mL is resistant.

¹¹Purity of this lot was assessed for 7 days at 37°C in an aerobic atmosphere on Tryptic Soy agar.





Date: 22 FEB 2017

Signature:

BEI Resources Authentication

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898