

NIH AIDS Reagent Program

20301 Century Boulevard Building 6, Suite 200 Germantown, MD 20874 USA

Phone: 240 686 4740 Fax: 301 515 4015 aidsreagent.org

DATA SHEET

Reagent:	$\dot{\mathbf{x}}$ SIV BK28 Infected HUT 78 Cells
Catalog Number:	173
Lot Number:	99067
Release Category:	D
Provided:	9 x 10 ⁶ cells/ml. Viability is 95%.
Propagation Medium:	RPMI 1640, 90%; fetal bovine serum, 10%.
Freeze Medium:	MEM with 50% fetal bovine serum, 90%; DMSO, 10%.
Growth Characteristics:	The cells grow in single cell suspension. Doubling time is 3 days; passage the cells every 3-4 days. Seeding 5×10^6 cells in a 75cm ² flask yields approximately $18-20 \times 10^6$ cells in a week. Addition of fresh, uninfected cells is not necessary to maintain the culture. When thawing, immediately place vial in a 37°C water bath, transfer cells to a sterile 12 x 45 mm centrifuge tube after 12-45 minutes, add fresh complete medium and centrifuge at 1200 RPM for 5 minutes. The cells do not grow in other media.
Sterility:	Negative for mycoplasma, bacteria and fungi.
Description:	HUT 78 Cells infected with SIV BK28.
Special Characteristics:	BK28 virus causes persistent infection and lymphadenopathy in rhesus macaques. Cellular atypia and giant cell formation occur in infected HUT 78 cells. The plasmid clone pBK28-SIV is available as Catalog #133.
Recommended Storage:	Liquid nitrogen.
Contributor:	Dr. James I. Mullins.

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.

References:	Kornfeld H, Riedel N, Viglianti GA, Hirsch V, Mullins JI. Cloning of HTLV-4 and its relation to simian and human immunodeficiency viruses. <i>Nature</i> 326 :610-613, 1987.
NOTE:	Acknowledgment for publications should read "The following reagent was obtained through the NIH AIDS Reagent Program, Division of AIDS, NIAID, NIH: SIV BK28 Infected HUT 78 Cells from Dr. James Mullins." Also include the reference cited above in any publications.

Last Updated July 02, 2018

ALL RECIPIENTS OF THIS MATERIAL MUST COMPLY WITH ALL APPLICABLE BIOLOGICAL, CHEMICAL, AND/OR RADIOCHEMICAL SAFETY STANDARDS INCLUDING SPECIAL PRACTICES, EQUIPMENT, FACILITIES, AND REGULATIONS. NOT FOR USE IN HUMANS.