

Certificate of Analysis for NR-52240

Acinetobacter baumannii, Strain MRSN 337038

Catalog No. NR-52240

This reagent is the tangible property of the U.S. Government.

Product Description:

Acinetobacter baumannii (A. baumannii), strain MRSN 337038 was isolated in 2015 from a human in Europe as part of a global surveillance program. NR-52240 was deposited as sensitive to colistin and tetracycline and resistant to amikacin, ceftazidime, ciprofloxacin, ceftriaxone, cefepime, gentamicin, imipenem, levofloxacin, meropenem, trimethoprim/sulfamethoxazole and ampicillin/sulbactam, and intermediately resistant to tobramycin. NR-52240 lot 70039376 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70039376 Manufacturing Date: 25SEP2020

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-negative rods	Gram-negative rods
Colony morphology	Report results	Circular, convex, entire, smooth and cream (Figure 1)
Growth at 44°C ± 2°C¹ 1 day in an aerobic atmosphere on Tryptic Soy agar	Growth	Growth
Motility Remel™ Motility Test Medium w/TTC Indicator for 1 day at 37°C in an aerobic atmosphere	Report results	Motile
VITEK® MS (MALDI-TOF)	A. baumannii	A. baumannii (99.9%)
Antibiotic Susceptibility Profile ^{2,3}		
Amikacin	Resistant	Resistant (> 256 µg/mL)
Ampicillin/sulbactam	Resistant	Resistant (32 µg/mL)
Cefepime	Resistant	Resistant (96 µg/mL)
Ceftriaxone	Resistant	Resistant (≥ 64 µg/mL)
Ceftazidime	Resistant	Resistant (≥ 64 µg/mL)
Ciprofloxacin	Resistant	Resistant (≥ 4 μg/mL)
Colistin	Sensitive	Sensitive (≤ 0.25 μg/mL)
Gentamicin	Resistant	Resistant (≥ 16 µg/mL)
Imipenem	Resistant	Resistant (> 8 μg/mL)
Levofloxacin	Resistant	Resistant (16 µg/mL)
Meropenem	Resistant	Resistant (> 8 µg/mL)
Trimethoprim/sulfamethoxazole	Resistant	Resistant (4 μg/mL)
Tobramycin	Intermediate	Sensitive (≤ 2 μg/mL) ⁴
Tetracycline	Sensitive	Resistant (> 128 µg/mL) ⁵
Genotypic Analysis		
Sequencing of 16S ribosomal RNA gene (~ 1480 base pairs)	≥ 99% sequence identity to A. baumannii, strain MRSN 337038 (GenBank: VHEY01000100.1)	100% sequence identity to A. baumannii, strain MRSN 337038 (GenBank: VHEY01000100.1)

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898



Certificate of Analysis for NR-52240

TEST	SPECIFICATIONS	RESULTS
Purity 7 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability	Growth	Growth

¹Growth at 44°C differentiates A. baumannii from A. calcoaceticus and A. pittii, which do not grow at 44°C.

⁵A. baumannii, strain MRSN 337038 was deposited as sensitive to tetracycline, but showed a MIC of > 128 μg/mL (interpreted as resistant) for tetracycline during QC testing. Testing was performed in duplicate.





/Heather Couch/ Heather Couch

27 MAY 2021

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

²Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

³Antibiotic susceptibility was tested using a combination of VITEK®2 GN82, Sensititre GNX2F AST and E-test strips.

⁴A. baumannii, strain MRSN 337038 was deposited as intermediate to tobramycin, but showed a MIC of ≤ 2 μg/mL (interpreted as sensitive) for tobramycin during QC testing. Testing was performed in duplicate.