

# **Certificate of Analysis for NR-52151**

### Acinetobacter baumannii, Strain MRSN 959

#### Catalog No. NR-52151

This reagent is the tangible property of the U.S. Government.

### **Product Description:**

Acinetobacter baumannii (A. baumannii), strain MRSN 959 was isolated in 2008 from a human wound in the United States as part of a global surveillance program. A. baumannii, strain MRSN 959 was deposited as sensitive to colistin and tetracycline and resistant to amikacin, ampicillin/sulbactam, cefepime, ceftazidime, ceftriaxone, ciprofloxacin, gentamicin, imipenem, levofloxacin, meropenem, tobramycin and trimethoprim/sulfamethoxazole. NR-52151 was produced by inoculation of BEI Resources seed lot 70038255 into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70048827 Manufacturing Date: 24NOV2021

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: Contact@BEIResources.org. We try to respond to feedback within 24 hours.

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-negative rods	Gram-negative rods
Colony morphology	Report results	Circular, convex, entire, smooth and cream
Growth at 44°C ± 2°C <sup>1</sup>	Growth	Growth
1 day in an aerobic atmosphere on Tryptic Soy agar		
Motility	Report results	Motile
BBL™ Motility Test Medium w/TTC Indicator for 1 day at 37°C in an aerobic atmosphere		
VITEK® MS (MALDI-TOF)	A. baumannii	A. baumannii (99.9%)
Antibiotic Susceptibility Profile <sup>2,3</sup>		
Amikacin	Resistant	Resistant (≥ 256 µg/mL)
Ampicillin/sulbactam	Resistant	Resistant (≥ 256 µg/mL)
Cefepime	Resistant	Resistant (≥ 256 µg/mL)
Ceftriaxone	Resistant	Resistant (≥ 64 µg/mL)
Ceftazidime	Resistant	Resistant (≥ 64 µg/mL)
Ciprofloxacin	Resistant	Resistant (≥ 4 µg/mL)
Colistin	Sensitive	Sensitive (≤ 0.25 μg/mL)
Gentamicin	Resistant	Resistant (≥ 16 µg/mL)
Imipenem	Resistant	Resistant (> 8 μg/mL)
Levofloxacin	Resistant	Resistant (12 µg/mL)
Meropenem	Resistant	Resistant (> 8 μg/mL)
Trimethoprim/sulfamethoxazole	Resistant	Resistant (> 4 μg/mL)
Tobramycin	Resistant	Resistant (≥ 256 µg/mL)
Tetracycline	Sensitive	Sensitive (3 µg/mL)
Genotypic Analysis		
Sequencing of 16S ribosomal RNA gene	≥ 99% sequence identity to	100% sequence identity to
(~ 1470 base pairs)	A. baumannii, strain MRSN 959 (GenBank: VHDS01000075.1)	A. baumannii, strain MRSN 959 (GenBank: VHDS01000075.1)
Purity	Growth consistent with expected	Growth consistent with expected
7 days at 37°C in an aerobic atmosphere with and without 5% CO <sub>2</sub> on Tryptic Soy agar	colony morphology	colony morphology

BEI Resources www.beiresources.org E-mail: <a href="mailto:contact@beiresources.org">contact@beiresources.org</a>
Tel: 800-359-7370

Fax: 703-365-2898



# **Certificate of Analysis for NR-52151**

TEST	SPECIFICATIONS	RESULTS
Viability	Growth	Growth

<sup>&</sup>lt;sup>1</sup>Growth at 44°C differentiates A. baumannii from A. calcoaceticus and A. pittii, which do not grow at 44°C.

/Sonia Bjorum Brower/ Sonia Bjorum Brower

11 OCT 2022

Technical Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Tel: 800-359-7370 Fax: 703-365-2898

<sup>&</sup>lt;sup>2</sup>Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

<sup>&</sup>lt;sup>3</sup>Antibiotic susceptibility was tested using a combination of VITEK® 2 GN82, Sensititre™ GNX2F AST and E-test strips.