

SUPPORTING INFECTIOUS DISEASE RESEARCH

Certificate of Analysis for NR-10322

Genomic DNA from Staphylococcus aureus, Strain HFH-30032

Catalog No. NR-10322

Product Description: Genomic DNA was obtained from a preparation of *Staphylococcus aureus* (*S. aureus*), strain HFH-30032 which was isolated from a human wound in Michigan, 2003.

Lot¹: 58859208 Manufacturing Date: 19NOV2009

TEST	SPECIFICATIONS	RESULTS
Sequencing of 16S Ribosomal RNA Gene (~ 1340 bp)	Consistent with S. aureus	Consistent with S. aureus
Agarose Gel Electrophoresis	High molecular weight chromosomal DNA	High molecular weight chromosomal DNA (Figure 1)
Concentration by PicoGreen® Measurement	4 to 6 μg in 25 to 100 μL per vial	5.4 μg in 52 μL per vial (104 μg/mL)
Functional Activity by PCR Amplification 16S ribosomal RNA gene	~ 1500 bp amplicon	~ 1500 bp amplicon
OD ₂₆₀ /OD ₂₈₀ Ratio	1.7 to 1.9	1.7
Bacterial Inactivation 10% of total yield plated on Tryptic Soy Agar with 5% sheep blood ^{2,3}	No viable bacteria detected	No viable bacteria detected

The bacterial preparation used for extraction of genomic DNA was produced by Brain Heart Infusion culture of the deposited material. After incubation for 24 hours at 37°C and aerobic atmosphere, genomic DNA was extracted using proprietary technology.

Date: 18 MAR 2010 **Signature:** Signature on File

Title: Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC® s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

Biodefense and Emerging Infections Research Resources Repository


www.beiresources.org

E-mail: contact@beiresources.org
Tel: 800-359-7370

²7 days at 37°C in an aerobic atmosphere.

³An extraction procedure was used that has been shown to consistently inactivate 100% of Gram-negative bacteria.

Certificate of Analysis for NR-10322



Lane 1: Invitrogen™ TrackIt 1 Kb Plus DNA Ladder™

Lane 2: 200 ng of NR-10322

NR-10322_58859208_18MAR2010

E-mail: contact@beiresources.org