SUPPORTING INFECTIOUS DISEASE RESEARCH

## Genomic RNA from Kilbourne F172: A/Leningrad/360/86 (HA, NA) x A/Puerto Rico/8/34 (H3N2), Reassortant X-91

## Catalog No. NR-10368

**Product Description:** Genomic RNA was isolated from a preparation of pooled allantoic fluid from specific-pathogen free embryonated chicken eggs infected with reassortant influenza A virus, A/Leningrad/360/86 (HA, NA) x A/Puerto Rico/8/34 (H3N2) (Kilbourne F172; X-91).

## Lot<sup>2</sup>: 58437301

## Manufacturing Date: 09JUL2009

| TEST   | SPECIFICATIONS     | RESULTS                           |
|--|--------------------|-----------------------------------|
| Sequencing of Influenza A Specific Region<br>Matrix gene (796 nucleotides)   | Influenza A virus  | Influenza A virus                 |
| Total RNA Content by RiboGreen <sup>®</sup> Measurement<br>(Viral, Cellular and Carrier)   | Report results     | 99 ng per 100 μL                  |
| Functional Activity by RT-PCR Amplification <sup>3,4</sup>   | ~ 1030 bp amplicon | ~ 1030 bp amplicon (See Figure 1) |
| Virus Inactivation<br>10% of total yield incubated in SPF embryonated<br>chicken eggs <sup>1</sup> and evaluated by hemagglutination<br>assay <sup>5</sup> | No virus detected  | No virus detected                 |
| Sodium Azide Content   | Report results     | 0.004%                            |

<sup>1</sup>10 to 11-day-old SPF Fertile Chicken Eggs acquired from B&E Eggs, York Springs, Pennsylvania

<sup>2</sup>Nucleic acid was extracted from a preparation of reassortant influenza A virus, A/Leningrad/360/86 (HA, NA) x A/Puerto Rico/8/34 (H3N2) (BEI Resources NR-3513; Lot 58253968) using a QIAamp<sup>®</sup> Viral RNA Mini kit (Qiagen 52906). <sup>3</sup>Amplified using 1:10 dilution of NR-10368 and One-Step RT-PCR Kit (Qiagen 210212)

<sup>4</sup>The primers are described in Hoffmann, E., et al. "Universal Primer Set for the Full-Length Amplification of All Influenza A Viruses." Arch. Virol. 146 (2001): 2275-2289. PubMed: 11811679.

<sup>5</sup>This extraction procedure has been shown to consistently inactivate 100% of influenza viruses.

Date: 26 JAN 2011

crothy C. Young Signature: 🔀

Title:


Technical Manager, BEI Authentication or designee

ATCC<sup>®</sup>, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC<sup>®</sup>'s knowledge.

ATCC<sup>®</sup> is a trademark of the American Type Culture Collection. You are authorized to use this product for research use only. It is not intended for human use.

800-359-7370 Fax: 703-365-2898 E-mail: contact@beiresources.org bieii resources

SUPPORTING INFECTIOUS DISEASE RESEARCH



Lane 1: 10-fold dilution of NR-10368 Lane 2: Invitrogen™ TrackIt™ 1 Kb Plus DNA Ladder

800-359-7370 Fax: 703-365-2898 E-mail: <u>contact@beiresources.org</u>