

## **Certificate of Analysis for NR-10082**

Genomic RNA from Kilbourne F110: A/equine/Prague/1/56 (HA) x A/Aichi/2/68 (NA) x A/Puerto Rico/8/34 (H7N2), Reassortant X-33, NA Deficient

Catalog No. NR-10082

**Product Description:** Genomic RNA was isolated from a preparation of pooled allantoic fluid from specific-pathogen free (SPF) embryonated chicken eggs<sup>1</sup> infected with reassortant influenza A virus, A/equine/Prague/1/56 (HA) x A/Aichi/2/68 (NA) x A/Puerto Rico/8/34 (H7N2) (Kilbourne F110; X-33, NA deficient).

Lot<sup>2</sup>: 58406990 Manufacturing Date: 29JUL2009

TEST	SPECIFICATIONS	RESULTS
Sequencing of Influenza A Specific Region Matrix gene (868 nucleotides)	Influenza A virus Identical to BEI Resources NR-3620	Influenza A virus Identical to BEI Resources NR-3620
Total RNA Content by RiboGreen® Measurement (Viral, Cellular and Carrier)	Report results	12.5 ng per 100 μL
Functional Activity by RT-PCR Amplification <sup>3,4</sup>	~ 1030 bp amplicon	~ 1030 bp amplicon (See Figure 1)
Virus Inactivation 10% of total yield incubated in SPF embryonated chicken eggs <sup>1</sup> and evaluated by hemagglutination assay <sup>5</sup>	No virus detected	No virus detected
Sodium Azide Content	Report results	0.0004%

<sup>10</sup> to 11-day-old SPF Fertile Chicken Eggs acquired from B&E Eggs, York Springs, Pennsylvania

**Date**: 25 JAN 2011 **Signature**:

Signature: Scrothy C. Young

Title: Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

Biodefense and Emerging Infections Research Resources Repository

Fax: 703-365-2898

E-mail: contact@beiresources.org

800-359-7370


<sup>&</sup>lt;sup>2</sup>Nucleic acid was extracted from a preparation of reassortant influenza A virus, Aequine/Prague/1/56 (HA) x A/Aichi/2/68 (NA) x A/Puerto Rico/8/34 (H7N2) (BEI Resources NR-3620; Lot 58253972) using a QIAamp<sup>®</sup> Viral RNA Mini kit (Qiagen 52906).

<sup>&</sup>lt;sup>3</sup>Amplified using 1:10 dilution of NR-10082 and One-Step RT-PCR Kit (Qiagen 210212)

<sup>&</sup>lt;sup>4</sup>The primers are described in Hoffmann, E., et al. "Universal Primer Set for the Full-Length Amplification of All Influenza A Viruses." <u>Arch. Virol.</u> 146 (2001): 2275-2289. PubMed: 11811679.

<sup>&</sup>lt;sup>5</sup>This extraction procedure has been shown to consistently inactivate 100% of influenza viruses.

## **Certificate of Analysis for NR-10082**



Lane 1: Invitrogen™ TrackIt™ 1 Kb DNA Ladder Lane 2: 10-fold dilution of NR-10082

Page 2 of 2

Fax: 703-365-2898

E-mail: contact@beiresources.org

800-359-7370