

Certificate of Analysis for NR-10326

Genomic DNA from Staphylococcus aureus, Strain HFH-30123

Catalog No. NR-10326

Product Description: Genomic DNA was obtained from a preparation of *Staphylococcus aureus* (*S. aureus*), strain HFH-30123 which was isolated from a human sinus in Michigan, 2004.

Lot¹: 59413901 Manufacturing Date: 25AUG2010

TEST	SPECIFICATIONS	RESULTS
Sequencing of 16S Ribosomal RNA Gene (~ 1470 base pairs)	Consistent with S. aureus	Consistent with S. aureus
Agarose Gel Electrophoresis	High molecular weight chromosomal DNA	High molecular weight chromosomal DNA (Figure 1)
Concentration by PicoGreen® Measurement	0.7 to 1.5 μg in 25 to 100 μL per vial	1.0 μg in 26 μL per vial (38 μg/mL)
Functional Activity by PCR Amplification 16S ribosomal RNA gene	~ 1500 bp amplicon	~ 1500 bp amplicon
OD ₂₆₀ /OD ₂₈₀ Ratio	1.7 to 1.9	1.8
Bacterial Inactivation 10% of total yield plated on Brain Heart Infusion Agar ^{2,3}	No viable bacteria detected	No viable bacteria detected

¹The bacterial preparation used for extraction of genomic DNA was produced by Brain Heart Infusion culture of the deposited material. After incubation for 24 hours at 37°C in an aerobic atmosphere, genomic DNA was extracted using proprietary technology.

Date: 11 MAY 2011 Signature:

Title: Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370


Fax: 703-365-2898

²7 days at 37°C in an aerobic atmosphere

³An extraction procedure was used that has been shown to consistently inactivate 100% of Gram-negative and Gram-positive bacteria.



Certificate of Analysis for NR-10326



Lane 1: Invitrogen™ TrackIt 1 Kb Plus DNA Ladder™

Lane 2: 200 ng of NR-10326

E-mail: contact@beiresources.org Tel: 800-359-7370

Fax: 703-365-2898