

## **Certificate of Analysis for HM-64D**

## Genomic DNA from Veillonella sp., Strain 3\_1\_44

Catalog No. HM-64D

**Product Description:** Genomic DNA was obtained from a preparation of *Veillonella* sp., strain 3 1 44.

Lot<sup>1,2</sup>: 60467432 Manufacturing Date: 06APR2012

| TEST   | SPECIFICATIONS  | RESULTS   |
|--|---|---|
| Sequencing of 16S Ribosomal RNA Gene<br>(~ 1110 base pairs)  | ≥ 99% identical to GenBank:<br>ADCV01000019 ( <i>Veillonella</i> sp.,<br>strain 3_1_44) | ≥ 99% identical to GenBank:<br>ADCV01000019 (Veillonella sp.,<br>strain 3_1_44) |
| Agarose Gel Electrophoresis  | High molecular weight chromosomal DNA   | High molecular weight chromosomal DNA (Figure 1)                                |
| Concentration by PicoGreen® Measurement  | 0.7 to 1.5 μg in 25 to 100 μL per vial  | 0.9 μg in 36 μL per vial (26 μg/mL)   |
| Functional Activity by PCR Amplification<br>16S ribosomal RNA gene   | ~ 1500 bp amplicon  | ~ 1500 bp amplicon  |
| OD <sub>260</sub> /OD <sub>280</sub> Ratio   | 1.7 to 2.0  | 1.8   |
| Bacterial Inactivation 10% of total yield inoculated on Tryptic Soy Agar with 5% defibrinated sheep blood <sup>3,4</sup> | No viable bacteria detected   | No viable bacteria detected   |

<sup>&</sup>lt;sup>1</sup>Quality control of HMP organisms used for DNA extraction is only performed to demonstrate that the material produced by BEI Resources is identical to the deposited material. It should not be considered a complete characterization of the deposited organism.

Date: 17 SEP 2012 Signature:

**Title:** Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898


<sup>&</sup>lt;sup>2</sup>The bacterial preparation used for extraction of genomic DNA was produced by culture of the deposited material. Genomic DNA was extracted using proprietary technology.

<sup>&</sup>lt;sup>3</sup>7 days at 37°C and anaerobic atmosphere

<sup>&</sup>lt;sup>4</sup>An extraction procedure was used that has been shown to consistently inactivate 100% of Gram-positive and Gram-negative bacteria.



## **Certificate of Analysis for HM-64D**



Lane 1: Invitrogen™ TrackIt 1 Kb Plus DNA Ladder™

Lane 2: 200 ng of HM-64D

E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898