

## **Certificate of Analysis for NR-28800**

## Salmonella enterica subsp. enterica, Strain SL486 (Serovar Heidelberg)

## Catalog No. NR-28800

Product Description: Salmonella enterica (S. enterica) subsp. enterica strain SL486 (also referred to as CVM35945 and N418), was isolated in 2005 from a human patient in Georgia, USA. Antibiotic resistance has not been reported for this strain.

Lot<sup>1</sup>: 61683857 Manufacturing Date: 19APR2013

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis	_	
Cellular morphology	Gram-negative rods	Gram-negative rods
Colony morphology <sup>2</sup>	Report results	Circular, low convex, entire, smooth and cream (Figure 1)
Analytical profile index (API <sup>®</sup> 20 E) <sup>3</sup>	Consistent with Salmonella spp.	Consistent with Salmonella spp.
Genotypic Analysis		
Sequencing of 16S ribosomal RNA gene (~ 960 base pairs)	Consistent with S. enterica	Consistent with <i>S. enterica</i> subsp. enterica
Riboprinter <sup>®</sup> Microbial Characterization System	Consistent with S. enterica	Consistent with S. enterica subsp. enterica serovar Heidelberg
Serogroup Analysis	Report results	Serogroup B <sup>4</sup>
Viability (post-freeze) <sup>2</sup>	Growth	Growth

<sup>&</sup>lt;sup>1</sup>NR-28800 was produced by inoculation of the deposited material into Nutrient broth and incubated 24 hours at 37°C in an aerobic atmosphere. Broth inoculum was added to kolles, which were incubated for 24 hours at 37°C in an aerobic atmosphere to produce this lot.

Figure 1



**Date: 14 JAN 2014** Signature:

> Title: Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

**BEI Resources** 

www.beiresources.org

E-mail: contact@beiresources.org Tel: 800-359-7370

Fax: 703-365-2898

<sup>&</sup>lt;sup>2</sup>24 hours at 37°C in an aerobic atmosphere on Nutrient agar

<sup>&</sup>lt;sup>3</sup>Additional serological tests are required to determine Salmonella species using API<sup>®</sup> 20 E.

<sup>&</sup>lt;sup>4</sup>Serogroup B contains serovar Heidelberg in addition to other serovars, see Ng, S. P., et al. "Detection and Serogroup Differentiation of Salmonella" spp. in Food within 30 Hours by Enrichment-Immunoassay with a T6 Monoclonal Antibody Capture Enzyme-Linked Immunosorbent Assay." Appl. Environ. Microbiol. 62 (1996): 2794-2302. PubMed: 8779567.