

Staphylococcal Enterotoxin B Toxoid, Chemically Inactivated from *Staphylococcus aureus* subsp. *aureus*

Catalog No. NR-4672

This reagent is the property of the U.S. government.

Product Description: Staphylococcal enterotoxin B (SEB) toxoid was obtained from a preparation of *Staphylococcus aureus* subsp. *aureus*, strain FDA 243 (ATCC[®] 14458[™]) and chemically inactivated with formaldehyde.

Lot: 57680184

Manufacturing Date: JUL2007

TEST	SPECIFICATIONS	RESULTS
Appearance	Clear with no particulate matter	Clear with no particulate matter
SDS-PAGE (SYPRO Orange densitometer scan)	Dominant band ~ 28 kDa > 95% pure	Dominant band of ~ 28 kDa > 95% pure (see Figure 1) ¹
SELDI-TOF Mass Spectrometry	Measured mass within 5% of expected mass based on amino acid sequence (~ 28,366 daltons)	Measured mass of untreated protein ~ 1% greater than expected mass Measured mass of treated protein ~ 4% greater than expected mass
SELDI-TOF Mass Spectrometry of Lys-C Proteolysis Products	> 50% of total residues accounted for in peptides of expected mass	57% of total residues accounted for in peptides of expected mass
Concentration by Bicinchoninic Acid Protein Assay	Report results	0.8 to 1.1 mg/mL
Functional Activity Western blot ² (see Figure 2) Staphylococcal enterotoxin B Staphylococcal enterotoxin B toxoid Carbonic anhydrase ELISA (see Figure 3) Staphylococcal enterotoxin B Staphylococcal enterotoxin B toxoid Negative control	Reactive Reactive Not reactive Report results Report results Not reactive	Reactive Reactive ¹ Not reactive 1:1,000,000 1:100,000 Not reactive
Cytotoxicity Human interferon-γ ELISPOT assay ^{3,4} Staphylococcal enterotoxin B Staphylococcal enterotoxin B toxoid Negative control (cells with no antigen)	Positive T-cell response Negative T-cell response Negative T-cell response	 > 300 spots per sample 20 spots per sample 22 spots per sample
Sterility	0.22 µm filter sterilized	0.22 µm filter sterilized
Absorbance Ratio (OD ₂₈₀ /OD ₂₆₀)	Report results	1.1
Endotoxin Content	Report results	12.5 EU/mg

¹The toxoid appears as a smear due to formaldehyde treatment of the protein.

²Using Anti-Staphylococcal Enterotoxin B antibody produced in rabbit (SIGMA[®] S9008).


³Using Human Interferon-γ ELISPOT Kit (Pierce ELHIFNG).

⁴This assay analyzes interferon-γ production at a single cell level. A spot forms at the site where secreted cytokine is bound; one spot is equal to one molecule-secreting cell.

Biodefense and Emerging Infections Research Resources Repository P.O. Box 4137 Manassas, VA 20108-4137 USA www.beiresources.org 800-359-7370 Fax: 703-365-2898 E-mail: <u>contact@beiresources.org</u>



Certificate of Analysis for NR-4672



Date: 14 DEC 2007

Signature: Signature on File

Title: Technical Manager, BEI Authentication or designee

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected by the contractor to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.



ATCC[®] is a trademark of the American Type Culture Collection. You are authorized to use this product for research use only. It is not intended for human use.

Biodefense and Emerging Infections Research Resources Repository P.O. Box 4137 Manassas, VA 20108-4137 USA www.beiresources.org 800-359-7370 Fax: 703-365-2898 E-mail: <u>contact@beiresources.org</u>