SUPPORTING INFECTIOUS DISEASE RESEARCH

#### Cryptococcus gattii, Strain CBS1930

#### Catalog No. NR-43209

**Product Description:** *Cryptococcus gattii* (*C. gattii*), strain CBS1930 was isolated from a goat in Aruba prior to the outbreak in Vancouver, British Columbia, Canada.

## Lot<sup>1,2</sup>: 61631742

## Manufacturing Date: 29MAR2013

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology <sup>3</sup>	Report results	Globose to ovoid, single or budding (Figure 1A)
Colony morphology <sup>3</sup>	Report results	Smooth, mucoid, entire and cream (Figure 1B)
Canavanine-glycine-bromthymol blue (CGB) differential medium <sup>4</sup>	Blue (C. gatti)	Blue ( <i>C. gatti</i> )
Genotypic Analysis		
Sequencing of partial 18S rRNA gene, internal transcribed spacer (ITS) 1, 5.8S rRNA gene,	Consistent with C. gattii	Consistent with <i>C. gattil⁵</i>
ITS 2, partial 28S rRNA (~ 540 base pairs) Sequencing of 26S rRNA gene (~ 620 base pairs)	Consistent with C. gattii	Consistent with <i>C. gattil<sup>5</sup></i>
Purity <sup>5</sup>		
Nutrient broth with 0.1% Yeast Extract at 25°C Nutrient broth with 0.1% Yeast Extract at 37°C	No bacterial growth No bacterial growth	No bacterial growth No bacterial growth
Viability (post-freeze) <sup>2</sup>	Growth	Growth

<sup>1</sup>NR-43209, lot 61631742, was produced by the depositor by incubation at 30°C in Yeast Peptone Dextrose medium overnight. The resultant growth was mixed with 30% glycerol to a final concentration of 15% and vialed.

<sup>2</sup>Quality control testing was performed on the vialed material by BEI Resources.


<sup>3</sup>4 days at 25°C in an aerobic atmosphere on Yeast Mold agar

<sup>4</sup>35 hours at 27°C in an aerobic atmosphere. CGB medium differentiates *C. gattii* from *C. neoformans* based on the ability of *C. gatti* isolates to grow in the presence of L-canavanine and to assimilate glycine as a sole carbon source, resulting in a blue color. *C. neoformans* isolates will show yellow to light-green on CGB medium. [McTaggart, L., et al. "Rapid Identification of *Cryptococcus neoformans* var. *grubii*, *C. neoformans* var. *neoformans*, and *C. gattii* by Use of Rapid Biochemical Tests, Differential Media, and DNA Sequencing." <u>J. Clin. Microbiol</u>. 2011 (49): 2522-2527. PubMed: 21593254.]

<sup>5</sup>Also consistent with *C. neoformans* 

<sup>6</sup>Clarity of broth was determined by visual inspection after 3 days at 25°C and 37°C in an aerobic atmosphere.





BEI Resources www.beiresources.org E-mail: <u>contact@beiresources.org</u> Tel: 800-359-7370 Fax: 703-365-2898 biei resources

# Certificate of Analysis for NR-43209

SUPPORTING INFECTIOUS DISEASE RESEARCH

Date: 28 JAN 2015

Signature: (

Jala.

Title:

Technical Manager, BEI Authentication or designee

ATCC<sup>®</sup>, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC<sup>®</sup>'s knowledge.

ATCC<sup>®</sup> is a trademark of the American Type Culture Collection. You are authorized to use this product for research use only. It is not intended for human use.

