

Certificate of Analysis for NR-48805

Genomic DNA from Neisseria meningitidis, Strain 9506

Catalog No. NR-48805

Product Description: Genomic DNA was obtained from a preparation of *Neisseria meningitidis* (*N. meningitidis*), strain 9506 which was isolated in 2004 from an infected patient in Brazil. *N. meningitidis*, strain 9506 was deposited as a member of serogroup C.

Lot¹: 62778964 Manufacturing Date: 11SEP2014

TEST	SPECIFICATIONS	RESULTS
Sequencing of 16S Ribosomal RNA Gene (~ 940 base pairs)	Consistent with N. meningitidis	Consistent with <i>N. meningitidis</i> ²
Agarose Gel Electrophoresis	High molecular weight chromosomal DNA	High molecular weight chromosomal DNA (Figure 1)
Concentration by PicoGreen® Measurement	0.7 to 1.5 μg in 25 to 100 μL per vial	1.0 μg in 27 μL per vial (37 μg/mL)
Functional Activity by PCR Amplification 16S ribosomal RNA gene	~ 1500 base pairs amplicon	~ 1500 base pairs amplicon
OD ₂₆₀ /OD ₂₈₀ Ratio	1.7 to 2.1	1.9
Bacterial Inactivation 10% of total yield plated on Chocolate agar ^{3,4}	No viable bacteria detected	No viable bacteria detected

¹The bacterial preparation used for extraction of genomic DNA was produced by broth culture of BEI Resources NR-32129 lot 62778969. Genomic DNA was extracted using proprietary technology.

Date: 08 APR 2015 Signature

Title: Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Tel: 800-359-7370 Fax: 703-365-2898


²100% identical *N. meningitidis*, strain 9506 (GenBank: ANRR01000003.1)

³7 days at 37°C in an aerobic atmosphere with 5% CO₂

⁴An extraction procedure was used that has been shown to consistently inactivate 100% of Gram-negative and Gram-positive bacteria.



Certificate of Analysis for NR-48805



Lane 1: Invitrogen™ TrackIt 1 Kb Plus DNA Ladder™ Lane 2: 200 ng of NR-48805