

Certificate of Analysis for NR-50053

Genomic DNA from Brucella abortus, Strain 870 (NCTC 10505)

Catalog No. NR-50053

Product Description: Genomic DNA was extracted from a preparation of *Brucella abortus* (*B. abortus*), strain 870 (NCTC 10505), biovar 6.

Lot¹: 63808663 Manufacturing Date: 30JAN2015

TEST	SPECIFICATIONS	RESULTS
Sequencing of 16S Ribosomal RNA Gene (~ 1380 base pairs)	Consistent with <i>B. abortus</i>	Consistent with <i>B. abortus</i> ²
Agarose Gel Electrophoresis	High molecular weight chromosomal DNA	High molecular weight chromosomal DNA (Figure 1)
Concentration by PicoGreen® Measurement	0.7 to 1.5 μg in 25 to 100 μL per vial	0.94 μg in 28 μL per vial (34 μg/mL)
Functional Activity by PCR Amplification 16S ribosomal RNA gene	~ 1500 base pair amplicon	~ 1500 base pair amplicon
OD ₂₆₀ /OD ₂₈₀ Ratio	1.7 to 2.1	1.6 ³
Bacterial Inactivation 10% of total yield plated on Tryptic Soy agar with 5% defibrinated sheep blood ^{4,5}	No viable bacteria detected	No viable bacteria detected

The bacterial preparation used for extraction of genomic DNA was produced by culture of BEI Resources NR-261 (lot 4741452). Genomic DNA was extracted using proprietary technology.

www.beiresources.org

Date: 17 FEB 2016 Signature:

BEI Authentication or designee

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

NR-50053 63808663 17FEB2016

²Also consistent with other *Brucella* species

³Lot 63808663 did not meet the OD₂₆₀/OD₂₈₀ ratio specifications, however, the integrity gel (Figure 1) shows intact genomic DNA and PCR amplification was successful. An OD₂₆₀/OD₂₈₀ ratio of 1.6 may indicate the presence of protein impurities.


⁴9 days at 37°C in an aerobic atmosphere with 5% CO₂

⁵An extraction procedure was used that has been shown to consistently inactivate 100% of Gram-negative and Gram-positive bacteria.



Certificate of Analysis for NR-50053

Figure 1: Agarose Gel Electrophoresis



Lane 1: Invitrogen™ Tracklt 1 Kb Plus DNA Ladder™ Lane 2: 200 ng of NR-50053

E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898