

## Genomic DNA from Escherichia coli, Strain RDEC-1

#### Catalog No. NR-3054

**Product Description:** Genomic DNA was isolated from a preparation of rabbit enteropathogenic *Escherichia coli*, strain RDEC-1, serotype O15:NM.

#### Lot<sup>1</sup>: 7642447

### Manufacturing Date: 22SEP2006

| TEST                                                                                      | SPECIFICATIONS                        | RESULTS                                          |
|-------------------------------------------------------------------------------------------|---------------------------------------|--------------------------------------------------|
| Sequencing of 16S Ribosomal RNA Gene (~ 475 bp)                                           | Consistent with Escherichia coli      | Consistent with Escherichia coli <sup>2</sup>    |
| Agarose Gel Electrophoresis                                                               | High molecular weight chromosomal DNA | High molecular weight chromosomal DNA (Figure 1) |
| Concentration by PicoGreen <sup>®</sup> Measurement                                       | 4 to 6 μg in 25 to 100 μL per vial    | 5.1 μg in 54 μL per vial (95 μg/mL)              |
| Functional Activity by PCR Amplification<br>16S ribosomal RNA gene                        | ~ 1500 bp amplicon                    | ~ 1500 bp amplicon                               |
| OD <sub>260</sub> /OD <sub>280</sub> Ratio                                                | 1.7 to 1.9                            | 1.8                                              |
| Bacterial Inactivation<br>10% of total yield plated on Trypticase Soy Agar <sup>3,4</sup> | No viable bacteria detected           | No viable bacteria detected                      |

<sup>1</sup>The bacterial preparation used for extraction of genomic DNA was produced by inoculation of BEI Resources NR-104 (Lot: 3681284) onto Trypticase Soy Agar. After incubation for 24 hours at 37°C, genomic DNA was extracted using proprietary technology.

<sup>2</sup>Also consistent with *Shigella* species.

<sup>3</sup>7 days at 37°C in an aerobic atmosphere.

<sup>4</sup>An extraction procedure was used that has been shown to consistently inactivate 100% of Gram-negative bacteria.

Date: 14 NOV 2007

# Signature: Signature on File


#### Title: Technical Manager, BEI Authentication or designee

ATCC<sup>®</sup>, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC<sup>®</sup>'s knowledge.

ATCC<sup>®</sup> is a trademark of the American Type Culture Collection. You are authorized to use this product for research use only. It is not intended for human use.

800-359-7370 Fax: 703-365-2898 E-mail: contact@beiresources.org





Lane 1: Invitrogen™ High DNA Mass™ Ladder Lane 2: 200 ng of NR-3054

800-359-7370 Fax: 703-365-2898 E-mail: contact@beiresources.org