

Escherichia coli, Strain H10407

Catalog No. NR-4

(Derived from ATCC® 35401™)

Product Description:

Escherichia coli (*E. coli*), strain H10407 was isolated from the feces of a patient with diarrhea in Bangladesh. NR-4 was produced by inoculation of BEI Resources seed lot 3561953 into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere with 5% CO₂ to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

Lot: 70058605

Manufacturing Date: 08FEB2023


TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology Motility (wet mount) VITEK® MS (MALDI-TOF)	Gram-negative rods Report results Report results <i>E. coli</i>	Gram-negative rods Circular, convex, entire, smooth and cream (Figure 1) Motile <i>E. coli</i> (99.9%)
Genotypic Analysis Next-Generation Sequencing (NGS) of Complete Genome ¹ Digital DNA-DNA hybridization (dDDH) ^{2,3}	≥ 99% sequence identity to <i>E. coli</i> , strain H10407 (GenBank: FN649414.1) ≥ 70% for species identification	100% sequence identity to <i>E. coli</i> , strain H10407 (GenBank: FN649414.1) <i>Shigella boydii</i> (85.4%) <i>Shigella sonnei</i> (85.2%) <i>Shigella flexneri</i> (84.1%) <i>Shigella dysenteriae</i> (81.2%) <i>E. coli</i> (74.7%)
Purity (post-freeze) 8 days at 37°C in an aerobic atmosphere with and without 5% CO ₂ on Tryptic Soy agar	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability (post-freeze)	Growth	Growth

¹The whole genome of *E. coli*, strain H10407 (~ 5.1 megabase pairs) was sequenced using the Illumina® MiSeq® system and was analyzed and assembled with Kraken version 2.1.2 and Unicycler v0.4.8.

²Relatedness between bacterial strains has traditionally been determined using DDH. For additional information, refer to Auch, A. F., et al. "Digital DNA-DNA Hybridization for Microbial Species Delineation by Means of Genome-to-Genome Sequence Comparison." *Stand. Genomic Sci.* 2 (2010): 117-134. PubMed: 21304684. dDDH analysis was performed using the Type (Strain) Genome Server.

³*Escherichia coli* and *Shigella* spp. are genetically related and cannot be differentiated by DDH. (Dif, G., et al. "In-Depth Genome-Based Analysis of *Shigella* spp. and *Escherichia* spp.: Resolving Ambiguities and Unveiling Phylogenetic Relationships." *Curr. Microbiol.* 82 (2025): 170. PubMed: 40045049.)

Figure 1: Colony Morphology



/Sonia Bjorum Brower/

Sonia Bjorum Brower

16 JUN 2025

Technical Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

