

**Leptospira kirschneri, Strain 200701401 (Serovar Bogvere)**

**Catalog No. NR-19942**

**Product Description:**

*Leptospira kirschneri* (*L. kirschneri*), strain 200701401 (Serovar Bogvere) was isolated in 2007 from a human in Guadeloupe, French West Indies. NR-19942 was produced by inoculation of BEI Resources seed lot 61906318 onto Ellinghausen-McCullough-Johnson-Harrison (EMJH) semisolid agar (0.15%), which was grown for 21 days at 30°C in an aerobic atmosphere. The material from the initial growth was passaged once in EMJH semisolid agar (0.15%) for 21 days at 30°C in an aerobic atmosphere to produce this lot. Quality control testing was completed under propagation conditions unless otherwise noted.

**Lot: 70077574**

**Manufacturing Date: 09SEP2025**

BEI Resources is committed to ensuring digital accessibility for people with disabilities. This Certificate of Analysis contains complex tables and may not be fully accessible. Please let us know if you encounter accessibility barriers and a fully accessible document will be provided: E-mail: [Contact@BEIResources.org](mailto:Contact@BEIResources.org). We try to respond to feedback within 24 hours.


TEST	SPECIFICATIONS	RESULTS
<b>Phenotypic Analysis</b> Cellular morphology Colony morphology  Motility (wet mount)	Spirochetes Growth below the soft agar surface (Dinger's disk) Motile	Spirochetes Growth below the soft agar surface (Dinger's disk) <sup>1</sup> (Figure 1) Motile
<b>Genotypic Analysis</b> Sequencing of 16S ribosomal RNA gene (~ 1450 base pairs)	Consistent with <i>L. kirschneri</i>	Consistent with <i>L. kirschneri</i> <sup>2</sup>
<b>Purity (post-freeze)</b> 21 days on Tryptic Soy agar with 5% defibrinated sheep blood at 37°C in an aerobic atmosphere with 5% CO <sub>2</sub>	No growth observed	No growth observed
<b>Viability (post-freeze)</b> Visual observation LIVE/DEAD® BacLight™ Bacterial Viability	Growth Green fluorescence visible	Growth <sup>3</sup> Green fluorescence visible <sup>3</sup> (Figure 2)

<sup>1</sup>Disk of dense growth below the soft agar surface (Dinger's disk) [Czekalowski, J. W., J. W. McLeod and J. Rodican. "The Growth and Respiration of *Leptospira* in Solid or Semi-Solid Media with Special Reference to Dinger's Phenomenon." *Br. J. Exp. Pathol.* 34 (1953): 588-595.] was evident after 21 days at 30°C in EMJH semisolid agar (0.15%).


<sup>2</sup>Also consistent with other *Leptospira* species

<sup>3</sup>Determined after 21 days incubation under cultivation conditions with LIVE/DEAD® BacLight™ Bacterial Viability Kit, 100x magnification (Invitrogen™ L34856). Cells with a compromised membrane that are dead or dying will stain red, while cells with an intact membrane will stain green.

**Figure 1: Dinger's Disk**



**Figure 2: LIVE/DEAD® BacLight™ Bacterial Viability**



/Sonia Bjorum Brower/

Sonia Bjorum Brower

17 FEB 2026

Technical Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

