

***Escherichia coli* Virulence Target EAF Primers**

Catalog No. NR-12204

This reagent is the tangible property of the U.S. Government.

For research use only. Not for human use.

Contributor:

NIH Biodefense and Emerging Infections Research Resources Repository, NIAID, NIH

Product Description:

Diarrheagenic *Escherichia coli* (*E. coli*) are classified into several pathogenic groups based on their virulence characteristics. NR-12204 contains forward and reverse primers that specifically amplify a region of the virulence target EAF, EPEC adherence factor, found on the plasmid pEAF of enteropathogenic *E. coli* (EPEC).¹

Material Provided:

Each vial contains approximately 80 µL of a mixture of forward and reverse primers in TE buffer (pH 7.0). The concentration is shown on the Certificate of Analysis.

Note: *E. coli* 12-Target Multiplex PCR 10X Buffer (BEI Resources NR-13440) will be provided with your shipment of NR-12204.

Packaging/Storage:

Primers were packaged aseptically in screw-capped plastic cryovials. The product is provided frozen on dry ice and should be stored at -60°C or colder upon arrival. Freeze-thaw cycles should be minimized.

Citation:

Acknowledgment for publications should read "The following reagent was obtained through the NIH Biodefense and Emerging Infections Research Resources Repository, NIAID, NIH: *Escherichia coli* Virulence Target EAF Primers, NR-12204."

Biosafety Level: 1

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the following publication: U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. [Biosafety in Microbiological and Biomedical Laboratories](#), 5th ed. Washington, DC: U.S. Government Printing Office, 2007; see www.cdc.gov/od/ohs/biosfty/bmb15/bmb15toc.htm.

Disclaimers:

You are authorized to use this product for research use only. It is not intended for human use.

Use of this product is subject to the terms and conditions of the BEI Resources Material Transfer Agreement (MTA). The MTA is available on our Web site at www.beiresources.org.

While BEI Resources uses reasonable efforts to include accurate and up-to-date information on this product sheet, neither ATCC® nor the U.S. Government make any warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. Neither ATCC® nor the U.S. Government warrants that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, use and disposal. ATCC® and the U.S. Government are not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, the U.S. Government, ATCC®, their suppliers and contributors to BEI Resources are not liable for damages arising from the misidentification or misrepresentation of products.

Use Restrictions:

This material is distributed for internal research, non-commercial purposes only. This material, its product or its derivatives may not be distributed to third parties. Except as performed under a U.S. Government contract, individuals contemplating commercial use of the material, its products or its derivatives must contact the contributor to determine if a license is required. U.S. Government contractors may need a license before first commercial sale. This material may be subject to third party patent rights.

References:

1. Kimata, K., et al. "Rapid Categorization of Pathogenic *Escherichia coli* by Multiplex PCR." [Microbiol. Immunol.](#) 49 (2005): 485-492. PubMed: 15965295.
2. Tobe, T., et al. "Complete DNA Sequence and Structural Analysis of the Enteropathogenic *Escherichia coli* Adherence Factor Plasmid." [Infect. Immun.](#) 67 (1999): 5455-5462. PubMed: 10496929.

ATCC® is a trademark of the American Type Culture Collection.



APPENDIX I

***E. coli* Virulence Target EAF Primers**

Recommended Reagents/Equipment

Reagent	Source	Catalog #
<i>E. coli</i> Virulence Target EAF Primers	BEI Resources	NR -12204
Positive Control Template, Genomic DNA from <i>E. coli</i> , Strain CDC (ATCC® 12807™)	BEI Resources	NR-3050
10X PCR Buffer	BEI Resources	NR-13440
GoTaq® Polymerase	Promega	M500B
dNTP Mix	Promega	U151
Molecular Biology Grade Water	ATCC®	60-2645

Reaction Mix¹

Reagent	Stock Concentration	Volume per Reaction (µL)
Molecular Biology Grade Water	---	19.2
10X PCR Buffer	10X	3
dNTP Mix	10 mM each	0.6
GoTaq® Polymerase	---	0.2
Primers ²	10 µM (each primer)	5
Template DNA	1 ng per µL	2
		Total – 30 µL

¹Reaction mix should be kept on bench-top cooler until ready for use.

²Primers are supplied at working stock concentrations.

Cycling Protocol

Cycle	# of Repeats	Step	Conditions
1	1	1	94°C for 5 minutes
2	30	1	94°C for 1 minute
		2	52°C for 1 minute
		3	72°C for 1 minute
3	1	1	72°C for 7 minutes
4	Indefinite	1	Hold at 4°C