

Certificate of Analysis for NR-4622

Genomic DNA from Salmonella enterica subsp. enterica, Strain LS1001

Catalog No. NR-4622

Product Description: Genomic DNA was obtained from a preparation of *Salmonella enterica* (*S. enterica*) subsp. *enterica*, strain LS1001, a derivative of strain 14028 which was isolated from chickens.

Lot¹: 58580366 Manufacturing Date: 02MAY2009

TEST	SPECIFICATIONS	RESULTS
Sequencing of 16S Ribosomal RNA Gene (~ 1410 bp)	Identical to BEI Resources NR-4341 Consistent with S. enterica	Identical to BEI Resources NR-4341 Consistent with S. enterica
Agarose Gel Electrophoresis	High molecular weight chromosomal DNA	High molecular weight chromosomal DNA (Figure 1)
Concentration by PicoGreen® Measurement	4 to 6 μg in 25 to 100 μL per vial	5.2 μg in 61 μL per vial (86 μg/mL)
Functional Activity by PCR Amplification 16S ribosomal RNA gene	~ 1500 bp amplicon	~ 1500 bp amplicon
OD ₂₆₀ /OD ₂₈₀ Ratio	1.7 to 1.9	1.8
Bacterial Inactivation 10% of total yield plated on Tryptic Soy Agar ^{2,3}	No viable bacteria detected	No viable bacteria detected

Salmonella enterica subsp. enterica, strain LS1001 is a derivative of strain 14028, which was isolated from chickens. The bacterial preparation used for extraction of genomic DNA was produced by Tryptic Soy Broth culture of the deposited material. After incubation for 24 hours at 37°C and aerobic atmosphere, genomic DNA was extracted using proprietary technology.

www.beiresources.org

Date: 01 OCT 2009 **Signature:** Signature on File

Title: Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

Fax: 703-365-2898


800-359-7370

²7 days at 37°C in an aerobic atmosphere.

³An extraction procedure was used that has been shown to consistently inactivate 100% of Gram-negative bacteria.

SUPPORTING INFECTIOUS DISEASE RESEARCH

Certificate of Analysis for NR-4622



Lane 1: Invitrogen™ TrackIt 1 Kb Plus DNA Ladder™

Lane 2: 200 ng of NR-4622

800-359-7370