

Certificate of Analysis for NR-2768

Genomic RNA from Dengue Virus Type 1, TH-Sman

Catalog No. NR-2768

Product Description: Genomic RNA was isolated from a preparation of cell lysate and supernatant from African green monkey (Vero) cells¹ infected with dengue virus type 1 (DEN-1), TH-Sman. RT-PCR amplification of a sequence of approximately 1200 nucleotides was successful with dilutions of up to 10-fold for this lot of NR-2768.

Lot²: 57835415 Manufacturing Date: 19AUG2009

TEST	SPECIFICATIONS	RESULTS
Sequencing of DEN-1 Specific Sequence (975 nucleotides)	Consistent with DEN-1	Consistent with DEN-1
Functional Activity by RT-PCR Amplification ³	~ 1200 bp amplicon	~ 1200 bp amplicon (see Figure 1)
Total RNA Content by RiboGreen [®] Measurement (Viral, Cellular, and Carrier)	Report results	456 ng per 100 μL
Virus Inactivation 10% of total yield inoculated on Vero cells and evaluated for immunofluorescence 1,4	No viable virus detected	No viable virus detected
Sodium Azide Content	Report results	0.004%

¹Vero cells: ATCC® CCL-81™

Date: 19 FEB 2010 **Signature:** Signature on File

Title: Technical Manager, BEI Authentication or designee

800-359-7370

Fax: 703-365-2898

NR-2768_57835415_19FEB2010

E-mail: contact@beiresources.org

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.


Biodefense and Emerging Infections Research Resources Repository

²Nucleic acid was extracted from a preparation of DEN-1, TH-Sman using a QIAGEN QIAamp[®] Viral RNA Mini Kit.

³Amplified using a 1:10 dilution of NR-2768 and a QIAGEN One-Step RT-PCR Kit (210212)

⁴This extraction procedure has been shown to consistently inactivate 100% of dengue viruses using a number of methods designed to detect the virus in cells: 1) cytopathic effect (100% plating, 10% plating, blind passage sub-culturing) and 2) indirect immunofluorescence of viral antigens.

Certificate of Analysis for NR-2768



Lane 1: Invitrogen™ TrackIt™ 1 Kb Plus DNA Ladder Lane 2: 10-fold dilution of NR-2768

Fax: 703-365-2898

E-mail: contact@beiresources.org

800-359-7370