

Certificate of Analysis for NR-13517

Genomic DNA from Acinetobacter baumannii, Isolate 7

Catalog No. NR-13517

Product Description: Genomic DNA was obtained from a preparation of *Acinetobacter baumannii* (*A. baumannii*) that was isolated from human nasotrach aspirate in 2008.

Lot¹: 58666820 Manufacturing Date: 09OCT2009

TEST	SPECIFICATIONS	RESULTS
Sequencing of 16S Ribosomal RNA Gene (~ 1460 bp)	Consistent with A. baumannii	Consistent with A. baumannii
Agarose Gel Electrophoresis	High molecular weight chromosomal DNA	High molecular weight chromosomal DNA (Figure 1)
Concentration by PicoGreen® Measurement	4 to 6 μg in 25 to 100 μL per vial	4.7 μg in 33 μL per vial (143 μg/mL)
Functional Activity by PCR Amplification 16S ribosomal RNA gene	~ 1500 bp amplicon	~ 1500 bp amplicon
OD ₂₆₀ /OD ₂₈₀ Ratio	1.7 to 1.9	1.7
Bacterial Inactivation 10% of total yield plated on Tryptic Soy Agar with 5% sheep blood ^{2,3}	No viable bacteria detected	No viable bacteria detected

The bacterial preparation used for extraction of genomic DNA was produced by Tryptic Soy Broth culture of the deposited material. After incubation for 24 hours at 37°C and aerobic atmosphere, genomic DNA was extracted using proprietary technology.

Date: 30 APR 2010 **Signature:** Signature on File

Title: Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

Biodefense and Emerging Infections Research Resources Repository

www.beiresources.org


Tel: 800-359-7370

Tel: 800-359-7370 Fax: 703-365-2898

²7 days at 37°C in an aerobic atmosphere

³An extraction procedure was used that has been shown to consistently inactivate 100% of Gram-negative bacteria

Certificate of Analysis for NR-13517



Lane 1: Invitrogen™ Tracklt 1 Kb Plus DNA Ladder™ Lane 2: 200 ng of NR-13517

NR-13517_58666820_30APR2010