

Certificate of Analysis for HM-162D

Genomic DNA from Staphylococcus aureus, Strain MN8

Catalog No. HM-162D

Product Description: Genomic DNA was obtained from a preparation of *Staphylococcus aureus* (*S. aureus*), strain MN8 was isolated from the vagina of a toxic shock syndrome patient.

Lot¹: 58954323 Manufacturing Date: 05JAN2010

TEST	SPECIFICATIONS	RESULTS
Sequencing of 16S Ribosomal RNA Gene (~ 1220 bp)	Identical to BEI Resources HM-162 Consistent with <i>S. aureus</i>	Identical to BEI Resources HM-162 Consistent with S. aureus
Agarose Gel Electrophoresis	High molecular weight chromosomal DNA	High molecular weight chromosomal DNA (Figure 1)
Concentration by PicoGreen® Measurement	0.7 to 1.5 μg in 25 to 100 μL per vial	1.1 μg in 37 μL per vial (29 μg/mL)
Functional Activity by PCR Amplification 16S ribosomal RNA gene	~ 1500 bp amplicon	~ 1500 bp amplicon
OD ₂₆₀ /OD ₂₈₀ Ratio	1.7 to 1.9	1.7
Bacterial Inactivation 10% of total yield plated on Tryptic Soy Agar with 5% defibrinated sheep blood ^{2,3}	No viable bacteria detected	No viable bacteria detected

¹The bacterial preparation used for extraction of genomic DNA was produced by Tryptic Soy Broth culture of the deposited material. After incubation for 24 hours at 37°C and aerobic atmosphere, genomic DNA was extracted using proprietary technology.

Date: 08 JUN 2010 **Signature:** Signature on file

Title: Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected by ATCC® to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

 $\mathsf{ATCC}^{\$}$ is a trademark of the American Type Culture Collection. You are authorized to use this product for research use only. It is not intended for human use.


E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

²7 days at 37°C in an aerobic atmosphere

³An extraction procedure was used that has been shown to consistently inactivate 100% of Gram-negative bacteria and Gram-positive bacteria.

Certificate of Analysis for HM-162D



Lane 1: Invitrogen™ TrackIt 1 Kb Plus DNA Ladder™

Lane 2: 200 ng of HM-162D

HM-162D_58954323_08JUN2010