

## Certificate of Analysis for HM-201D

## Genomic DNA from Enterococcus faecalis, Strain TX0104

Catalog No. HM-201D

**Product Description:** Genomic DNA was obtained from a preparation of *Enterococcus faecalis* (*E. faecalis*), strain TX0104.

Lot<sup>1</sup>: 59079887 Manufacturing Date: 18MAR2010

TEST	SPECIFICATIONS	RESULTS
Sequencing of 16S Ribosomal RNA Gene (~ 1410 base pairs)	Consistent with E. faecalis	Consistent with E. faecalis <sup>2</sup>
Agarose Gel Electrophoresis	High molecular weight chromosomal DNA	High molecular weight chromosomal DNA (Figure 1)
Concentration by PicoGreen® Measurement	0.75 to 1.25 μg in 25 to 100 μL per vial	1.2 μg in 28 μL per vial (43 μg/mL)
Functional Activity by PCR Amplification 16S ribosomal RNA gene	~ 1500 bp amplicon	~ 1500 bp amplicon
OD <sub>260</sub> /OD <sub>280</sub> Ratio	1.7 to 1.9	1.9
Bacterial Inactivation 10% of total yield plated on Tryptic Soy Agar <sup>3,4</sup>	No viable bacteria detected	No viable bacteria detected

<sup>&</sup>lt;sup>1</sup>The bacterial preparation used for extraction of genomic DNA was produced by Tryptic Soy Broth culture of the deposited material. After incubation for 24 hours at 37°C and aerobic atmosphere, genomic DNA was extracted using proprietary technology.

Date: 10 OCT 2010 Signature:

Title: Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected by ATCC® to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

Biodefense and Emerging Infections Research Resources Repository

www.beiresources.org

Tel: 800-359-7370


HM-201D\_59079887\_10OCT2010

<sup>&</sup>lt;sup>2</sup>Also consistent with *Enterococcus lactis* 

<sup>&</sup>lt;sup>3</sup>7 days at 37°C in an aerobic atmosphere

<sup>&</sup>lt;sup>4</sup>An extraction procedure was used that has been shown to consistently inactivate 100% of Gram-positive and Gram-negative bacteria.

## **Certificate of Analysis for HM-201D**



Lane 1: Invitrogen™ TrackIt 1 Kb Plus DNA Ladder™

Lane 2: 200 ng of HM-201D

HM-201D\_59079887\_10OCT2010

E-mail: contact@beiresources.org