

## **Certificate of Analysis for NR-13410**

## pART Companion Vector for Gene Expression in *Mycobacterium tuberculosis* and *Escherichia coli*

Catalog No. NR-13410

This reagent is the tangible property of the U.S. Government.

**Product Description:** pART is a companion vector that can be used in either *Escherichia coli* or *Mycobacterium tuberculosis*. The pART vector contains origins of replication for both organisms, Mycobacteriophage L5 integrase CDS and *att*P site, a *Mycobacterium smegmatis* acetamidase promoter region and the *Streptomyces coelicolor* tetracycline repressor CDS.

Lot: 59310230 Manufacturing Date: 16JUN2010

QC testing was performed by Colorado State University under the TB Vaccine Testing and Research Materials Contract (NIH). The Colorado State University documentation is attached.

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected by the contractor to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC<sup>®</sup> is a trademark of the American Type Culture Collection. You are authorized to use this product for research use only. It is not intended for human use.

NR-13410\_59310230\_15OCT2010

E-mail: contact@beiresources.org

## Recombinant Plasmid Quality Control Record

Plasmid designation	pART
BEI Product Item Nu	mber <u>NR-13410</u>
BEI Lot Number	59310230
CSU Lot Number	10.pART.6.16
Notebook/Pgs	Noebook #5/Page 23 (NKG)
Notebook detail	
Media used	<u>LB</u>
Culture size	250 mL
Growth conditions:	Temp 37 Time 18 hrs Shaker speed 200
Plasmid prep type (mi (Cat. No. 12643)	ni/maxi, kit name or protocol) Qiagen HiSpeed Plasmid Midi Kit
Plasmid prep detail:	Midi prep Qia100 tip and elution conditions
Strain used to produce	e plasmid TOP10
E. coli ori? Y/N	<u>Y</u>
Contains Mycobacteri	al ori? Y/N Y
Final concentration	158.97 ng/μL
Total Stocks119_	<del>-</del>
Buffer TE	-
Method used for quan	tifying nanodrop
QC gel – N/A	
Restriction enzymes u	sed in QC analysis N/A
Expected size of restri	ction fragments
	•

Created on 7/9/2010 11:16:00 AM Page 2

Vector	<u>N/A</u>		
Insert	<u>N/A</u>		
Other	N/A		
Gel description	file number, % agarose, buffer	N/A	
Molecular weight markers N/A			
Generated by Nutrolan May		Date 8/19/10	
Supervisor	e ice,s	Date 8 (20710	

Form 4.2.09.KMD