

Certificate of Analysis for NR-10079

Genomic RNA from Kilbourne F3: A/duck/Singapore/645/97 (H5N3) Mutant, High (Hy) Yield

Catalog No. NR-10079

Product Description: Genomic RNA was isolated from a preparation of pooled allantoic fluid from specific-pathogen free (SPF) embryonated chicken eggs¹ infected with a high (hy) yield mutant (Kilbourne F3) of influenza A virus, A/duck/Singapore/645/97 (H5N3).

Lot²: 58776474 Manufacturing Date: 08SEP2009

TEST	SPECIFICATIONS	RESULTS
Sequencing of Influenza A Specific Region Matrix gene (896 nucleotides)	Influenza A virus Identical to BEI Resources NR-3599	Influenza A virus Identical to BEI Resources NR-3599
Total RNA Content by RiboGreen [®] Measurement (Viral, Cellular and Carrier)	Report results	22 ng per 100 μL
Functional Activity by RT-PCR Amplification ^{3,4}	~ 1030 bp amplicon	~ 1030 bp amplicon (See Figure 1)
Virus Inactivation 10% of total yield incubated in SPF embryonated chicken eggs ¹ and evaluated by hemagglutination assay ⁵	No virus detected	No virus detected
Sodium Azide Content	Report results	0.0004%

¹⁰ to 11-day-old SPF Fertile Chicken Eggs acquired from B&E Eggs, York Springs, Pennsylvania

Date: 24 JAN 2011

Signature: Coung C. Young

Title: Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

Biodefense and Emerging Infections Research Resources Repository

Fax: 703-365-2898

E-mail: contact@beiresources.org

800-359-7370


²Nucleic acid was extracted from a preparation of high yield mutant of influenza A virus, A/duck/Singapore/645/97 (H5N3) (BEI Resources NR-3599; Lot 58253970) using a QIAamp® Viral RNA Mini kit (Qiagen 52906).

Amplified using 1:10 dilution of NR-10079 and One-Step RT-PCR Kit (Qiagen 210212)

⁴The primers are described in Hoffmann, E., et al. "Universal Primer Set for the Full-Length Amplification of All Influenza A Viruses." <u>Arch. Virol.</u> 146 (2001): 2275-2289. PubMed: 11811679.

⁵This extraction procedure has been shown to consistently inactivate 100% of influenza viruses.

Certificate of Analysis for NR-10079



Lane 1: 10-fold dilution of NR-10079

Lane 2: Invitrogen™ TrackIt™ 1 Kb DNA Ladder

Fax: 703-365-2898

E-mail: contact@beiresources.org

800-359-7370