

Certificate of Analysis for HM-121D

Genomic DNA from Streptococcus salivarius, Strain SK126

Catalog No. HM-121D

Product Description: Genomic DNA was obtained from a preparation of *Streptococcus salivarius* (*S. salivarius*), strain SK126.

Lot^{1,2}: 59980221 Manufacturing Date: 10JUN2011

TEST	SPECIFICATIONS	RESULTS
Sequencing of 16S Ribosomal RNA Gene (~ 1420 base pairs)	≥ 99% identical to depositor's sequence Consistent with <i>S. salivarius</i>	Pending Consistent with S. salivarius ³
Agarose Gel Electrophoresis	High molecular weight chromosomal DNA	High molecular weight chromosomal DNA (Figure 1)
Concentration by PicoGreen® Measurement	0.7 to 1.5 μg in 25 to 100 μL per vial	1.1 μg in 62 μL per vial (17 μg/mL)
Functional Activity by PCR Amplification 16S ribosomal RNA gene	~ 1500 bp amplicon	~ 1500 bp amplicon
OD ₂₆₀ /OD ₂₈₀ Ratio	1.7 to 2.0	1.8
Bacterial Inactivation 8% of total yield plated on Tryptic Soy Agar with 5% defibrinated sheep blood ^{4,5}	No viable bacteria detected	No viable bacteria detected

Quality control of HMP organisms used for DNA extraction is only performed to demonstrate that the material produced by BEI Resources is identical to the deposited material. It should not be considered a complete characterization of the deposited organism.

Date: 15 NOV 2011 **Signature:**

Title: Technical Manager, BEI Authentication or designee

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected by ATCC[®] to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources
www.beiresources.org
Tel: 800-359-7370

HM-121D_59980221_15NOV2011

²The bacterial preparation used for extraction of genomic DNA was produced by Tryptic Soy Broth culture of the deposited material. After incubation for 24 hours at 37°C in an aerobic atmosphere with 5% CO₂, genomic DNA was extracted using proprietary technology.


³Also consistent with *Streptococcus thermophilus*

⁴7 days at 37°C and 5% CO₂

⁵An extraction procedure was used that has been shown to consistently inactivate 100% of Gram-positive and Gram-negative bacteria.



Certificate of Analysis for HM-121D



Lane 1: Invitrogen™ TrackIt 1 Kb Plus DNA Ladder™

Lane 2: 200 ng of HM-121D

HM-121D_59980221_15NOV2011