b|**e**|**i** resources

SUPPORTING INFECTIOUS DISEASE RESEARCH

Schistosoma mansoni, Microsatellite CA11-1 Reverse Primer

Catalog No. NR-29395

This reagent is the tangible property of the U.S. Government.

For research use only. Not for human use.

Contributor:

Matty Knight, Principal Investigator, Biomedical Research Institute, Rockville, MD (NIH-NIAID Contract HHSN272201000005I)

Manufacturer:

Eurofins MWG Operon

Product Description:

NR-29395 is a twenty nucleotide reverse primer designed to amplify the microsatellite CA11-1 region from *Schistosoma mansoni* (*S. mansoni*) when paired with the CA11-1 forward primer (NR-29394). The sequence of the CA11-1-R 20-mer is 5'-CGATTAAAATGTACGTTTCC-3'.^{1,2} Please see Appendix I for general PCR procedure details.

Material Provided:

Each vial contains approximately 30 μ L of forward primer in nuclease-free distilled water. The concentration is shown on the Certificate of Analysis.

Packaging/Storage:

Primers were packaged aseptically in screw-capped plastic cryovials. The product is provided frozen on dry ice and should be stored at -20°C or colder upon arrival. Freeze-thaw cycles should be minimized. Note: For long-term storage it is strongly recommended that primers are kept in TE buffer (10 mM Tris, 1 mM EDTA), pH 8.0 at -20°C or colder.

Citation:

Acknowledgment for publications should read "The following reagent was obtained through BEI Resources, NIAID, NIH: *Schistosoma mansoni,* Microsatellite CA11-1 Reverse Primer, NR-29395."

Biosafety Level: 1

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the following publication: U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. <u>Biosafety in</u> <u>Microbiological and Biomedical Laboratories</u>. 5th ed. Washington, DC: U.S. Government Printing Office, 2009; see <u>www.cdc.gov/biosafety/publications/bmbl5/index.htm</u>.

Disclaimers:

You are authorized to use this product for research use only. It is not intended for human use.

Use of this product is subject to the terms and conditions of the BEI Resources Material Transfer Agreement (MTA). The MTA is available on our Web site at <u>www.beiresources.org</u>.

While BEI Resources uses reasonable efforts to include accurate and up-to-date information on this product sheet, neither ATCC[®] nor the U.S. Government makes any warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. Neither ATCC[®] nor the U.S. Government warrants that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, use and disposal. $ATCC^{\circledast}$ and the U.S. Government are not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, the U.S. Government, $ATCC^{\circledast}$, their suppliers and contributors to BEI Resources are not liable for damages arising from the misidentification or misrepresentation of products.

Use Restrictions:

This material is distributed for internal research, noncommercial purposes only. This material, its product or its derivatives may not be distributed to third parties. Except as performed under a U.S. Government contract, individuals contemplating commercial use of the material, its products or its derivatives must contact the contributor to determine if a license is required. U.S. Government contractors may need a license before first commercial sale. This material may be subject to third party patent rights.

References:

- Gower, C. M., et al. "Development and Application of an Ethically and Epidemiologically Advantageous Assay for the Multi-Locus Microsatellite Analysis of *Schistosoma mansoni.*" <u>Parasitology</u> 134 (2007): 523-536. PubMed: 17096873.
- Blair, L., Webster, J. P., and G. C. Barker "Isolation and Characterization of Polymorphic Microsatellite Markers in *Schistosoma mansoni* from Africa." <u>Mol. Ecol Notes</u> 1 (2001): 93-95.

ATCC[®] is a trademark of the American Type Culture Collection.



BEI Resources www.beiresources.org E-mail: <u>contact@beiresources.org</u> Tel: 800-359-7370 Fax: 703-365-2898 **b**|**e**|**i** resources

SUPPORTING INFECTIOUS DISEASE RESEARCH

APPENDIX I

S. mansoni Microsatellite Primers

Reagent	Source	Catalog #		
S. mansoni microsatellite primers (forward and reverse)	BEI Resources	NR-29379 to NR-29395		
Genomic DNA from <i>S. mansoni</i> ¹	BEI Resources	NR-28910 to NR-28912		
10X PCR Buffer	No Manufacturer Recommended	N/A		
Taq [®] Polymerase	No Manufacturer Recommended	N/A		
dNTP Mix	No Manufacturer Recommended	N/A		
Molecular Biology Grade Water	No Manufacturer Recommended	N/A		

commanded Passants/Equipment

¹Primers can also be used with other *S. mansoni* nucleic acids.

Reaction Mix¹

Reagent	Stock Concentration	Volume per Reaction (µL)
Molecular Biology Grade Water		16.5
10X PCR Buffer	10X	2.5
dNTP Mix	5 mM each	1
Taq [®] Polymerase	5 Units per μL	1
Forward and Reverse Primers ²	10 µM (each primer)	1
Template DNA	25 ng per µL	2
		Total – 25 μL

¹Reaction mix should be kept on bench-top cooler until ready for use.

²Primers are supplied at a concentration of 100 µM and need to be diluted to the working stock concentrations.

Cycle	# of Repeats	Step	Conditions
1	1	1	94°C for 15 minutes
2	10	1	94°C for 30 seconds
		2	68°C for 90 seconds (decrease temperature 1°C per cycle)
3	20	1	94°C for 30 seconds
		2	58°C for 90 seconds
4	1	1	72°C for 60 seconds
5	1	1	60°C for 30 minutes