

## **Certificate of Analysis for HM-259D**

## Genomic DNA from Treponema denticola, Strain F0402

Catalog No. HM-259D

**Product Description:** Genomic DNA was obtained from a preparation of *Treponema denticola* (*T. denticola*), strain F0402.

Lot<sup>1,2</sup>: 59980913 Manufacturing Date: 14SEP2012

TEST	SPECIFICATIONS	RESULTS
Sequencing of 16S Ribosomal RNA Gene (~ 1430 base pairs)	≥ 99% identical to GenBank: ADEC01000002 ( <i>T. denticola</i> , strain F0402)	≥ 99% identical to GenBank: ADEC01000002 ( <i>T. denticola</i> , strain F0402)
Agarose Gel Electrophoresis	High molecular weight chromosomal DNA	High molecular weight chromosomal DNA (Figure 1)
Concentration by PicoGreen® Measurement	0.7 to 1.5 μg in 25 to 100 μL per vial	1.0 μg in 44 μL per vial (24 μg/mL)
Functional Activity by PCR Amplification 16S ribosomal RNA gene	~ 1500 bp amplicon	~ 1500 bp amplicon
OD <sub>260</sub> /OD <sub>280</sub> Ratio	1.7 to 2.0	1.9
Bacterial Inactivation 10% of total yield inoculated in Modified New Oral Spirochete (NOS) Broth <sup>3-5</sup>	No viable bacteria detected	No viable bacteria detected

<sup>&</sup>lt;sup>1</sup>Quality control of HMP organisms used for DNA extraction is only performed to demonstrate that the material produced by BEI Resources is identical to the deposited material. It should not be considered a complete characterization of the deposited organism.

Date: 21 MAY 2013 Signature:

**Title:** Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

Page 1 of 2

ATCC® is a trademark of the American Type Culture Collection.

**BEI Resources** 

www.beiresources.org

You are authorized to use this product for research use only. It is not intended for human use.

E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

<sup>&</sup>lt;sup>2</sup>The bacterial preparation used for extraction of genomic DNA was produced by culture of the deposited material. Genomic DNA was extracted using proprietary technology.


<sup>&</sup>lt;sup>3</sup>7 days at 37°C in an anaerobic atmosphere (85% N<sub>2</sub>:5% CO<sub>2</sub>:10% H<sub>2</sub>)

<sup>&</sup>lt;sup>4</sup>Inoculated broth was analyzed for bacterial growth microscopically

<sup>&</sup>lt;sup>5</sup>An extraction procedure was used that has been shown to consistently inactivate 100% of Gram-positive and Gram-negative bacteria.



## **Certificate of Analysis for HM-259D**



Lane 1: 200 ng of HM-259D Lane 2: Invitrogen™ TrackIt 1 Kb Plus DNA Ladder™

E-mail: contact@beiresources.org Tel: 800-359-7370