

## **Certificate of Analysis for NR-28797**

## Salmonella enterica subsp. enterica, Strain SL483 (H7606) (Serovar Agona)

Catalog No. NR-28797

**Product Description**: Salmonella enterica (S. enterica) subsp. enterica strain SL483 (also referred to as strain H7606), was isolated in 2003 from a human patient in Wisconsin, USA. This strain is reported to be antibiotic-susceptible.

Lot<sup>1</sup>: 61683854 Manufacturing Date: 26APR2013

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Report results	Gram-negative rod
Colony morphology <sup>2</sup>	Report results	Circular, convex, entire, smooth and cream (Figure 1)
Analytical profile index (API <sup>®</sup> 20 E) <sup>3</sup>	Consistent with Salmonella spp.	Consistent with Salmonella spp.
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 770 base pairs)	Consistent with S. enterica	Consistent with S. enterica subsp.
Riboprinter® Microbial Characterization System	Consistent with S. enterica	Consistent with S. enterica subsp. enterica serovar Agona
Serogroup Analysis	Report results	Serogroup B <sup>5</sup>
Viability (post-freeze) <sup>2</sup>	Growth	Growth

NR-28797 was produced by inoculation of the deposited material into Nutrient broth and incubated 24 hours at 37°C in an aerobic atmosphere. Broth inoculum was added to kolles, which were incubated for 24 hours at 37°C in an aerobic atmosphere to produce this lot.

Figure 1



**Date:** 06 NOV 2013

Signature:

Title: Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC® s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

**BEI Resources** 

www.beiresources.org

E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

<sup>&</sup>lt;sup>2</sup>24 hours at 37°C in an aerobic atmosphere on Nutrient agar

<sup>&</sup>lt;sup>3</sup>Additional serological tests are required to determine Salmonella species using API<sup>®</sup> 20 E.

<sup>&</sup>lt;sup>4</sup>≥ 99% identical to GenBank: CP001138 (S. enterica subsp. enterica, serovar Agona, strain SL483)

<sup>&</sup>lt;sup>5</sup>Serogroup B contains serovar Agona in addition to other serovars [Ng, S. P., et al. "Detection and Serogroup Differentiation of *Salmonella* spp. in Food within 30 Hours by Enrichment-Immunoassay with a T6 Monoclonal Antibody Capture Enzyme-Linked Immunosorbent Assay." <u>Appl. Environ. Microbiol.</u> 62 (1996): 2794-2302. PubMed: 8779567.].