

Certificate of Analysis for NR-43785

N1 Neuraminidase (NA) Protein with N-Terminal Histidine Tag from Influenza Virus, A/Brisbane/59/2007 (H1N1), Recombinant from Baculovirus

Catalog No. NR-43785

This reagent is the tangible property of the U.S. Government.

Product Description: A recombinant form of the N1 neuraminidase (NA) protein from influenza A virus A/Brisbane/59/2007 (H1N1) containing an N-terminal histidine tag was produced in Sf9 insect cells using a baculovirus expression vector system and purified by nickel affinity chromatography.

Lot: 62190122 Manufacturing Date: 20NOV2013

TEST	SPECIFICATIONS	RESULTS
Appearance	Clear and colorless	Clear and colorless
Purity by SDS-PAGE Densitometry Scan	Protein band of interest represents ≥ 95% of total staining intensity	Protein band of ~ 55 kDa accounts for ~ 95% of total staining intensity (Figure 1)
Identification by Western Blot Analysis Polyclonal anti-N1 NA ¹ Monoclonal anti-histidine tag ²	Reactive Reactive	Reactive (Figure 2A) Reactive (Figure 2B)
Concentration by Bradford Assay ³	Report results	922 µg per mL
Final Product Quantity per vial Volume per vial	Report results Report results	101 μg 110 μL
Functional Activity Neuraminidase activity in fluorescent enzymatic assay	Report results	4.5 × 10 ⁷ relative fluorescence units per hour per mg protein ⁴
Filtration	0.22 µm filtered	0.22 µm filtered

¹BEI Resources NR-19260, Ferret Hyperimmune Sera to Influenza A/Brisbane/59/2007 (H1N1) (1:1000 dilution)

Date: 18 FEB 2014

Signature: Milharl Q. Complex

Title:

Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898

²Clontech (Cat. No. 631212) (1:1000 dilution)


³Using BSA as a standard

⁴Using serial dilutions of NR-43785 and 0.15 mM 2'-(4-methylumbelliferyl)-α-D-N-acetylneuraminic acid (4-MUNANA), Sigma (Cat. No. M8639), as described in Wetherall, N.T., et al. "Evaluation of Neuraminidase Enzyme Assays Using Different Substrates to Measure Susceptibility of Influenza Virus Clinical Isolates to Neuraminidase Inhibitors: Report of the Neuraminidase Inhibitor Susceptibility Network." J. Clin. Microbiol. 41 (2003): 742-750. PubMed: 12574276.



Certificate of Analysis for NR-43785

Figure 1: SDS-PAGE




Lane 1: SeeBlue® Plus2 Pre-Stained standard

Lane 2: NR-43785, 1 μg Lane 2: NR-43785, 2 μg

Figure 2: Western Blot Analysis


A. Polyclonal Anti-N2 NA



Lane 1: MagicMark™ XP Protein Standard

Lane 2: NR-43785, 0.5 µg

B. Monoclonal Anti-Histidine Tag



Lane 1: MagicMark™ XP Protein Standard

Lane 2: NR-43785, 0.5 μg