SUPPORTING INFECTIOUS DISEASE RESEARCH

Complementary DNA from Adult Schistosoma japonicum, Chinese Strain

Catalog No. NR-48637

This reagent is the tangible property of the U.S. Government.

Product Description: Complementary DNA (cDNA) was synthesized from RNA extracted from adult *Schistosoma japonicum*, Chinese strain.

Lot^{1,2}: 62506675

Manufacturing Date: 19FEB2014

TEST	SPECIFICATIONS	RESULTS
Concentration	Report results	1 µg in 20 µL per vial (0.05 µg/µL)
OD ₂₆₀ /OD ₂₈₀ Ratio	1.70 to 2.00	1.74
Qualification by RT-PCR Amplification of 28S ribosomal RNA gene ³	~ 290 base pair amplicon	~ 290 base pair amplicon (Figure 1)

¹QC testing was performed by the Biomedical Research Institute, Rockville, MD (NIH-NIAID Contract HHSN272201000005))

²cDNA was synthesized from total RNA by ProtoScript II[®] First Strand cDNA Synthesis Kit (New England BioLabs[®], Inc.), according to the manufacturer's instructions.

³Primers were designed to amplify the nucleotide region 39 to 326 of *Schistosoma mansoni (S. mansoni)* 28S ribosomal RNA gene (GenBank: Z46503.1). Cross-amplification of the 28S gene from other *Schistosoma* species has been observed.

Figure 1: Amplification of 28S Ribosomal RNA Gene for Qualification by RT-PCR bp Lane 1: 100 base pair ladder Lane 2: 290 base pair amplicon from S. mansoni 28S ribosomal RNA gene 500 400 300 200 100 ~ 290 base pairs

Date: 13 MAY 2014

Signature:

Title: Technical Manager, BEI Authentication or designee

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected by the contributor to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.



ATCC[®] is a trademark of the American Type Culture Collection. You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: <u>contact@beiresources.org</u> Tel: 800-359-7370 Fax: 703-365-2898