b|**e**|**i** resources

SUPPORTING INFECTIOUS DISEASE RESEARCH

Borrelia burgdorferi, Signature-Tagged Mutagenesis Library Clone T06TC166 (Gene BB_0431)

Catalog No. NR-23736

Product Description: Borrelia burgdorferi (B. burgdorferi), strain B31 5A18NP1 STM library clone T06TC166 was produced by signature-tagged mutagenesis (STM) of the BB_0431 gene.

Lot¹: 62782865

Manufacturing Date: 29JUL2014

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology ²	Spirochete	Spirochete
Motility (wet mount)	Motile	Motile
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1420 base pairs)	Consistent with <i>B. burgdorferi</i>	Consistent with <i>B. burgdorferi</i> ³
Confirmation of BBE02 Disruption ²	Growth in the presence of kanamycin	Growth in the presence of kanamycin
Confirmation of STM ²	Growth in the presence of gentamicin	Growth in the presence of gentamicin
Purity (post-freeze) ⁴	No growth observed	No growth observed
Viability (post-freeze)		
Visual observation	Growth	Growth ²
LIVE/DEAD [®] <i>Bac</i> Light [™] Bacterial Viability ⁵	Green fluorescence visible	Green fluorescence visible (Figure 1) ⁵


NR-23736 was produced by inoculation of the deposited material into Revised Barbour-Stoenner-Kelly medium supplemented with 200 µg/mL kanamycin and 40 µg/mL gentamicin and grown 20 days at 32°C in an aerobic atmosphere. The material from the initial growth was passaged once in Barbour-Stoenner-Kelly medium supplemented with 200 µg/mL kanamycin and 40 µg/mL gentamicin for 20 days at 32°C in an aerobic atmosphere to produce this lot.

²6 days at 32°C in an aerobic atmosphere Revised Barbour-Stoenner-Kelly medium supplemented with 200 μg/mL kanamycin and 40 μg/mL gentamicin

³>99.9% identical to GenBank: AE000783 (*B. burgdorferi,* strain B31)

⁴7 days at 37°C and 5% CO₂ on Tryptic Soy agar with 5% defibrinated sheep blood

⁵Determined after 6 days incubation under cultivation conditions with LIVE/DEAD[®] BacLight[™] Bacterial Viability Kit, 100x magnification (Invitrogen[™] L34856). Cells with a compromised membrane that are dead or dying will stain red, while cells with an intact membrane will stain green.



Date: 26 SEP 2014

Signature: Aald and

Title:

Technical Manager, BEI Authentication or designee

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

 $ATCC^{\otimes}$ is a trademark of the American Type Culture Collection. You are authorized to use this product for research use only. It is not intended for human use.



BEI Resources www.beiresources.org E-mail: contact@beiresources.org Tel: 800-359-7370 Fax: 703-365-2898