

## **Certificate of Analysis for NR-46889**

Staphylococcus aureus subsp. aureus, Strain JE2, Transposon Mutant NE346 (SAUSA300\_1346)

Catalog No. NR-46889


**Product Description:** Staphylococcus aureus (S. aureus) subsp. aureus, transposon mutant NE346 was derived from S. aureus subsp. aureus, strain JE2. Mutagenesis occurred through the use of the mariner-based transposon bursa aurealis resulting in an erythromycin-resistant deletion strain of JE2. S. aureus subsp. aureus, transposon mutant NE346 was created by disruption of SAUSA300\_1346, which encodes for a DnaQ family exonuclease/DinG family helicase that resembles the epsilon subunit of DNA polymerase III in Escherichia coli. Strain JE2 is a plasmid-cured derivative of strain LAC that was isolated in 2002 from a skin and soft tissue infection of an inmate in the Los Angeles County Jail in California, USA.

Lot<sup>1</sup>: 62959109 Manufacturing Date: 19SEP2014

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology <sup>2</sup> Motility (wet mount)	Gram-positive cocci Report results Report results	Gram-positive cocci Circular, convex, entire, smooth and cream (Figure 1) Non-motile
Confirmation of Transposon Insertion <sup>3</sup>	Resistant to erythromycin	Resistant to erythromycin
Viability (post-freeze) <sup>2</sup>	Growth	Growth

<sup>&</sup>lt;sup>1</sup>NR-46889 was produced by inoculation of the deposited material into Tryptic Soy broth with 5 μg/mL erythromycin and incubated for 21 hours at 37°C in an aerobic atmosphere.

Figure 1: Colony Morphology



BEI Resources
www.beiresources.org

E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

<sup>&</sup>lt;sup>2</sup>24 hours at 37°C in an aerobic atmosphere on Tryptic Soy agar with 5 µg/mL erythromycin

<sup>&</sup>lt;sup>3</sup>Prior to initiating work, it is recommended that the presence and location of the transposon is confirmed. Gene specific primers should be paired with either the "Upstream" primer (5'-CTCGATTCTATTAACAAGGG-3') for transposons in the "plus" orientation or the "Buster" primer (5'-GCTTTTTCTAAATGTTTTTTAAGTAAATCAAGTAC-3') for transposons in the "minus" orientation. For additional information refer to Fey, P. D., et al. "A Genetic Resource for Rapid and Comprehensive Phenotype Screening of Nonessential *Staphylococcus aureus* Genes." MBio 4 (2013): e00537-12. PubMed: 23404398.



## **Certificate of Analysis for NR-46889**

Date: 21 OCT 2014

Signature: (

**BEI** Resources Authentication

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898