

Certificate of Analysis for NR-46544

Staphylococcus aureus subsp. aureus, Strain JE2, Transposon Mutant NE1 (SAUSA300_1327)

Catalog No. NR-46544


Product Description: *Staphylococcus aureus* (*S. aureus*) subsp. *aureus*, transposon mutant NE1 was derived from *S. aureus* subsp. *aureus*, strain JE2. Mutagenesis occurred through the use of the *mariner*-based transposon *bursa aurealis* resulting in an erythromycin resistant deletion strain of JE2. *S. aureus* subsp. *aureus*, transposon mutant NE1 was created by disruption of SAUSA300_1327, which encodes for a cell surface protein. Strain JE2 is a plasmid-cured derivative of strain LAC that was isolated in 2002 from a skin and soft tissue infection of an inmate in the Los Angeles County Jail in California, USA.

Lot¹: 63025336 Manufacturing Date: 14OCT2014

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology ² Motility (wet mount)	Gram-positive cocci Report results Non-motile	Gram-positive cocci Circular, convex, entire, smooth and cream (Figure 1) Non-motile
Confirmation of Transposon Insertion ³	Resistant to erythromycin	Resistant to erythromycin
Purity (post-freeze) ⁴	Consistent with S. aureus	Consistent with S. aureus
Viability (post-freeze) ²	Growth	Growth

¹NR-46544 was produced by inoculation of the deposited material into Tryptic Soy broth with 5 μg/mL erythromycin and incubated for 22 hours at 37°C in an aerobic atmosphere.

 $^{^4}$ Purity of this lot was assessed for 7 days at 37°C in an aerobic atmosphere on Tryptic Soy agar with 5% defibrinated sheep blood.



BEI Resources
www.beiresources.org

E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

²24 hours at 37°C in an aerobic atmosphere on Tryptic Soy agar with 5 μg/mL erythromycin

³Prior to initiating work, it is recommended that the presence and location of the transposon is confirmed. Gene specific primers should be paired with either the "Upstream" primer (5'-CTCGATTCTATTAACAAGGG-3') for transposons in the "plus" orientation or the "Buster" primer (5'-GCTTTTTCTAAATGTTTTTTAAGTAAATCAAGTAC-3') for transposons in the "minus" orientation. For additional information, refer to Fey, P. D., et al. "A Genetic Resource for Rapid and Comprehensive Phenotype Screening of Nonessential *Staphylococcus aureus* Genes." MBio 4 (2013): e00537-12. PubMed: 23404398.



Certificate of Analysis for NR-46544

Date: 05 DEC 2014

Signature: (

Title:

Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources
www.beiresources.org

E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

NR-46544 63025336 05DEC2014