

## **Certificate of Analysis for NR-41916**

## Klebsiella pneumoniae, Strain BIDMC 2A

## Catalog No. NR-41916

Product Description: Klebsiella pneumoniae (K. pneumoniae), strain BIDMC 2A was isolated in 2010 from a blood culture in Boston, Massachusetts, USA.

Lot<sup>1</sup>: 62907100 Manufacturing Date: 11SEP2014

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology Colony morphologies <sup>2,3</sup>	Gram-negative bacilli	Gram-negative bacilli
Colony morphologies <sup>2,3</sup>	Report results	Colony type 1: Circular, convex, entire, smooth and cream (Figure 1)
		Colony type 2: Circular, low convex, entire, smooth, translucent and cream (Figure 1)
Motility (wet mount)	Report results	Non-motile
VITEK <sup>®</sup> MS (MALDI-TOF)	Consistent with K. pneumoniae	Consistent with K. pneumoniae
Genotypic Analysis		
Sequencing of 16S ribosomal RNA gene (~ 1480 base pairs) <sup>3</sup>	Consistent with K. pneumoniae	Consistent with K. pneumoniae
Riboprinter® Microbial Characterization System	Consistent with K. pneumoniae	Consistent with K. pneumoniae
Purity (post-freeze) <sup>4</sup>	Consistent with K. pneumoniae	Consistent with K. pneumoniae
Viability (post-freeze) <sup>2</sup>	Growth	Growth

<sup>&</sup>lt;sup>1</sup>K. pneumoniae, strain BIDMC 2A was deposited by Andrew B. Onderdonk, Director, Clinical Microbiology, Pathology, Brigham and Women's Hospital, Boston, Massachusetts, USA. NR-41916 was produced by inoculation of the deposited material into Nutrient broth and incubated for 22 hours at 37°C in an aerobic atmosphere. Broth inoculum was added to Nutrient agar kolles which were grown 24 hours at 37°C in an aerobic atmosphere to produce this lot.

Figure 1



Colony type 1 Colony type 2

**BEI Resources** www.beiresources.org E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

<sup>&</sup>lt;sup>2</sup>27 hours at 37°C in an aerobic atmosphere on Nutrient agar.

<sup>&</sup>lt;sup>3</sup>Two colony types were observed. Plating of the individual colony types showed that they did not revert to a single colony type. The 16S ribosomal RNA gene of each colony type was sequenced and found to be consistent with the other colony type and with K. pneumoniae.

<sup>&</sup>lt;sup>4</sup>Purity of this lot was assessed for 7 days at 37°C in an aerobic atmosphere on Nutrient agar.



## **Certificate of Analysis for NR-41916**

**Date:** 11 DEC 2014

Signature: (

Title:

Technical Manager, BEI Authentication or designee

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898