

Certificate of Analysis for NR-19133

Streptococcus pneumoniae, Strain GA19998

Catalog No. NR-19133

Product Description: *Streptococcus pneumoniae* (*S. pneumoniae*), strain GA19998 was isolated in 2004 from the blood of a patient with pneumonia in Georgia, USA. *S. pneumoniae*, strain GA19998 was deposited as a member of serotype 19A.

Lot¹: 62743341 Manufacturing Date: 02JUL2014

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-positive cocci	Gram-positive cocci
Colony morphologies (Figure 1) ^{2,3}	Report results	Colony type 1: Circular,
colony marphologica (rigulary)		umbilicate, undulate, opaque
		and gray
		Colony type 2: Circular, low
		convex, entire, opaque,
		smooth and gray
Hemolysis on blood agar ²	α-hemolytic	α-hemolytic
Motility (wet mount)	Report results	Non-motile
Catalase	Negative	Negative
VITEK [®] 2 (GP Card)	Consistent with S. pneumoniae	Consistent with S. pneumoniae
VITEK [®] MS (MALDI-TOF)	Consistent with S. pneumoniae	Consistent with S. pneumoniae
Antibiotic Susceptibility Profile ⁴		
Benzylpenicillin ⁵	Intermediate	Intermediate (4 µg/mL)
Ciprofloxacin ⁶	Report results	Resistant (3.0 µg/mL)
Levofloxacin ⁵	Sensitive	Sensitive (1 µg/mL)
Gatifloxacin ⁶	Report results	Sensitive (0.25 µg/mL)
Ofloxacin ⁵	Report results	Sensitive (2 µg/mL)
Erythromycin ⁵	Resistant	Resistant (≥ 1 µg/mL)
Telithromycin ⁵	Sensitive	Sensitive (0.5 µg/mL)
Clindamycin ⁶	Sensitive	Sensitive (0.25 µg/mL)
Quinupristin/dalfopristin ⁶	Sensitive	Sensitive (0.38 µg/mL)
Linezolid ⁵	Sensitive	Sensitive (≤ 2 μg/mL)
Vancomycin ⁵	Sensitive	Sensitive (≤ 1 µg/mL)
Tetracycline ⁵	Sensitive	Sensitive (≤ 1 µg/mL)
Rifampicin ⁶	Report results	Sensitive (0.064 µg/mL)
Trimethoprim/sulfamethoxazole ⁵	Resistant	Resistant (≥ 320 µg/mL)
Amoxicillin ⁵	Intermediate	Intermediate (4 µg/mL)
Meropenem ⁵	Intermediate	Intermediate (0.5 µg/mL)
Cefotaxime ⁵	Resistant	Resistant (≥ 4 µg/mL)
Cefuroxine ⁶	Resistant	Resistant (32.0 µg/mL)
Cefaclor ⁶	Report results	Resistant (256 µg/mL)
Ceftriaxone ⁵ Chloramphenicol ⁵	Resistant Sensitive	Resistant (≥ 4 µg/mL) Sensitive (≤ 2 µg/mL)
•	Sensitive	Sensitive (> 2 µg/mL)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene	Consistent with S. pneumoniae	Consistent with S. pneumoniae ⁷
(~ 1480 base pairs)	Consistent with S. prieumonide	Consistent with S. pheumoniae
Riboprinter® Microbial Characterization System	Consistent with S. pneumoniae	Consistent with S. pneumoniae
Purity (post-freeze) ⁸	Consistent with S. pneumoniae	Consistent with S. pneumoniae
Viability (post-freeze) ²	Growth	Growth

BEI Resources www.beiresources.org E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898



Certificate of Analysis for NR-19133

¹S. pneumoniae, strain GA19998 (also referred to as SPAR61) was deposited by Scott T. Chancey, Ph.D., Division of Infectious Diseases, Department of Medicine, Emory University, Atlanta, Georgia, USA. NR-19133 was produced by inoculation of the deposited material into Tryptic Soy broth and incubated for 24 hours at 37°C in an aerobic atmosphere with 5% CO₂. Broth inoculum was added to Tryptic Soy agar with 5% defibrinated sheep blood kolles which were grown 24 hours at 37°C in an aerobic atmosphere with 5% CO₂ to produce this lot.

²21 hours at 37°C in an aerobic atmosphere with 5% CO₂ on Tryptic Soy agar with 5% defibrinated sheep blood

³Two colony types were observed. Plating of the individual colony types showed that colony type 1 reverted back to a mixed colony type whereas colony type 2 remained true. The 16S ribosomal RNA gene of each colony type was sequenced and found to be consistent with *S. pneumoniae*.


⁴Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S22 (2012)

⁵Tested by VITEK[®] AST-GP74 card (VITEK[®] 2, version 5.04)

⁶Tested by bioMérieux E-test[®]: cefuroxime (catalog number 506958), ciprofloxacin (412310), clindamycin (412314), cefaclor (504550), gatifloxacin (530250), rifampicin (412449) and quinupristin/dalfopristin (528750)

⁷≥ 99.8% identical to S. pneumoniae, strain GA19998 (GenBank: ALCS01000005.1)

⁸Purity of this lot was assessed for 7 days on Tryptic Soy agar with 5% defibrinated sheep blood at 37°C in an aerobic atmosphere with 5% CO₂.



Date: 29 DEC 2014

Signature(

Title:

Technical Manager, BEI Authentication or designee

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898