SUPPORTING INFECTIOUS DISEASE RESEARCH

Cryptococcus gattii, Strain AIR265a

Catalog No. NR-43220

Product Description: Cryptococcus gattii (C. gattii), strain AIR265a is the progeny of a genotypic cross between C. gattii strains R265 and Alg166 and is one strain of a congenic pair.

Lot^{1,2}: 61631753

Manufacturing Date: 29MAR2013

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology ³	Report results	Globose to ovoid, single or budding
Colony morphology ³	Report results	(Figure 1A) Smooth, mucoid, entire and cream (Figure 1B)
Canavanine-glycine-bromthymol blue (CGB) differential medium ⁴	Blue (<i>C. gatti</i>)	Blue (<i>C. gatti</i>)
Genotypic Analysis Sequencing of partial 18S rRNA gene, internal transcribed spacer (ITS) 1, 5.8S rRNA gene, ITS 2, partial 28S rRNA (~ 510 base pairs)	Consistent with C. gattii	Consistent with <i>C. gattii⁵</i>
Sequencing of 26S rRNA (~ 510 base pairs) (~ 620 base pairs)	Consistent with C. gattii	Consistent with <i>C. gattii</i> ⁵
Purity ⁶ Nutrient broth with 0.1% Yeast Extract at 25°C Nutrient broth with 0.1% Yeast Extract at 37°C	No bacterial growth No bacterial growth	No bacterial growth No bacterial growth
Viability (post-freeze) ²	Growth	Growth

¹NR-43220, lot 61631753, was produced by the depositor by incubation at 30°C in Yeast Peptone Dextrose medium overnight. The resultant growth was mixed with 30% glycerol to a final concentration of 15% and vialed.

²Quality control testing was performed at BEI Resources.


³2 days at 25°C in an aerobic atmosphere on Yeast Mold agar

⁴35 hours at 27°C in an aerobic atmosphere. CGB medium differentiates *C. gattii* from *C. neoformans* based on the ability of *C. gatti* isolates to grow in the presence of L-canavanine and to assimilate glycine as a sole carbon source, resulting in a blue color. *C. neoformans* isolates will show yellow to light-green on CGB medium. [McTaggart, L., et al. "Rapid Identification of *Cryptococcus neoformans* var. *grubii*, *C. neoformans* var. *neoformans*, and *C. gattii* by Use of Rapid Biochemical Tests, Differential Media, and DNA Sequencing." J. Clin. Microbiol. 2011 (49): 2522-2527. PubMed: 21593254.]

⁵Also consistent with *C. neoformans*

⁶Clarity of broth was determined by visual inspection after 2 days at 25°C and 37°C in an aerobic atmosphere.

Figure 1: Cellular morphology (A) and colony morphology (B)



biei resources

Certificate of Analysis for NR-43220

SUPPORTING INFECTIOUS DISEASE RESEARCH

Date: 24 FEB 2015

Signature:

Title: Technical Manager, BEI Authentication or designee

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection. You are authorized to use this product for research use only. It is not intended for human use.

