

Escherichia coli K-12, Strain IM08B

Catalog No. NR-49806

Product Description: *Escherichia coli* (*E. coli*) K-12, strain IM08B contains the *hsdM* (methylase) and *hsdS* (specificity) genes from *Staphylococcus aureus* NRS384 clonal complex 8 (CC8). This insertion mutant was produced in *E. coli* K-12, strain DC10B (Δdcm). *E. coli* K-12, strain IM08B was deposited as resistant to streptomycin.

Lot¹: 63849802

Manufacturing Date: 18NOV2015

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology Colony morphology ² Motility (wet mount) VITEK [®] MS (MALDI-TOF)	Gram-negative rods Report results Report results Consistent with <i>E. coli</i>	Gram-negative rods Circular, low convex, entire, smooth and gray (Figure 1) Motile <i>E. coli</i> (99.9%)
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1480 base pairs) Riboprinter [®] Microbial Characterization System	≥ 99% sequence identity to <i>E. coli</i> K-12 strain (GenBank: NZ_CP014225.1) ≥ 85% <i>E. coli</i>	99.5% sequence identity to <i>E. coli</i> K-12 strain (GenBank: NZ_CP014225.1) 97% <i>E. coli</i>
Analysis of <i>hsdMS</i> by PCR Assay³ <i>hsdM2</i> and <i>hsdS2</i> (CC8) <i>hsdS1</i> (CC8)	~ 3400 base pair amplicon ~ 1770 base pair amplicon	~ 3400 base pair amplicon ~ 1770 base pair amplicon
Analysis of <i>hsdMS</i> by Sequence Analysis³ <i>hsdM2</i> and <i>hsdS2</i> (~ 1710 base pairs) <i>hsdS1</i> (~ 1390 base pairs)	Consistent with depositor sequence Consistent with depositor sequence	Consistent with depositor sequence Consistent with depositor sequence
Confirmation of Streptomycin Resistance²	Growth	Growth
Purity (post-freeze)⁴	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability (post-freeze)²	Growth	Growth

¹NR-49806 was produced by inoculation of the deposited material into Tryptic Soy broth with 25 µg/mL streptomycin and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar with 25 µg/mL streptomycin kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot.

²1 day at 37°C in an aerobic atmosphere on Tryptic Soy agar with 25 µg/mL streptomycin

³PCR primers used for amplification were IM434 forward primer 5'-ACTTTCTTTAAGGCTTAGAGTCAAGC-3', IM435 reverse primer 5'-TTTAACGCCACGTTCACTCTTTGC-3', 179 forward primer 5'-CGCCATTTATACAGGAAAAGCCTA-3' and 180 reverse primer 5'-GTTACCTTCTCTATAGAGAGTGGTG-3'. For additional information, refer to Monk, I., et al. "Complete Bypass of Restriction Systems for Major *Staphylococcus aureus* Lineages." *MBio*. 26 (2015): e00308-15. PubMed: 26015493.

⁴Purity of this lot was assessed for 7 days at 37°C in an aerobic atmosphere on Tryptic Soy agar with 5% defibrinated sheep blood.

Figure 1: Colony Morphology



Date: 15 APR 2016

Signature: 

BEI Resources Authentication

ATCC[®], on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC[®]'s knowledge.

ATCC[®] is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

