

Certificate of Analysis for NR-50400

Enterobacter cloacae complex, Strain BEI10

Catalog No. NR-50400

Product Description: Enterobacter cloacae complex (E. cloacae complex), strain BEI10 is from an unknown origin.

Lot¹: 64391840 Manufacturing Date: 30MAR2016

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-negative rods	Gram-negative rods
Colony morphology ²	Report results	Circular, convex, entire, smooth and
	·	cream (Figure 1)
Motility (wet mount)	Report results	Motile
Beta-lactamase ³	Report results	Positive
VITEK® 2 Compact (GN card)	≥ 90% probability of being	E. cloacae complex
,	E. cloacae complex	(98% probability) ⁴
Antibiotic Susceptibility Profile		
VITEK® (AST-GN84 Card)5,6		
Amoxicillin/Clavulanic Acid	Report results	Resistant (≥ 32 µg/mL)
Piperacillin/Tazobactam	Report results	Resistant (≥ 128 µg/mL)
Cefazolin	Report results	Resistant (≥ 64 µg/mL)
Ceftriaxone	Report results	Resistant (≥ 64 µg/mL)
Cefepime	Report results	Resistant (≥ 64 µg/mL)
Aztreonam	Report results	Resistant (≥ 64 µg/mL)
Ertapenem	Report results	Resistant (≥ 8 μg/mL)
Imipenem	Report results	Resistant (≥ 16 µg/mL)
Meropenem	Report results	Resistant (≥ 16 µg/mL)
Gentamicin	Report results	Resistant (≥ 16 µg/mL)
Ciprofloxacin	Report results	Resistant (≥ 4 µg/mL)
Levofloxacin	Report results	Resistant (≥ 8 μg/mL)
Tetracycline	Report results	Resistant (≥ 16 µg/mL)
Nitrofurantoin	Report results	Resistant (≥ 512 µg/mL)
Trimethoprim/Sulfamethoxazole	Report results	Resistant (≥ 320 µg/mL)
Etest® antibiotic test strips ⁷		
Ampicillin ⁸	Report results	Resistant (≥ 256 µg/mL)
Genotypic Analysis		
Sequencing of 16S ribosomal RNA gene	≥ 99% sequence identity to	100% sequence identity to
(~ 870 base pairs)	E. cloacae complex type strain	E. cloacae complex type strain
	(Genbank: NR_118568.1)	(Genbank: NR_118568.1) ⁹
Purity (post-freeze) ¹⁰	Consistent with expected colony	Consistent with expected colony
	morphology	morphology
Viability (post-freeze) ²	Growth	Growth
		<u> </u>

¹NR-50400 was produced by inoculation of the deposited material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot.

BEI Resources www.beiresources.org

E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

²1 day at 37°C in an aerobic atmosphere on Tryptic Soy agar

³The production of beta-lactamase was detected using a Cefinase™ Paper Disc (BBL™ 231650).

⁴Percent probabilities above 90% indicate a close match to the typical biochemical pattern for the given organism, with a percent probability of 99% being a perfect match between the test reaction pattern and the unique biochemical pattern of the given organism or organism group. For additional information, please refer to O'Hara, C.M. and J. M. Miller. "Evaluation of the VITEK 2 ID-GNB Assay for Identification of Members of the Family Enterobacteriaceae and Other Nonenteric Gram-Negative Bacilli and Comparison with the VITEK GNI+ Card." J. Clin. Microbiol. 41 (2003): 2096-2101. PubMed: 12734254.

⁵Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S22 (2012)

⁶No results were obtained for ampicillin, ampicillin

⁷¹ day at 37°C in an aerobic atmosphere on Mueller Hinton agar



Certificate of Analysis for NR-50400

⁸For ampicillin (bioMérieux Etest[®] 412252), a MIC ≤ 8 μg/mL is sensitive, a MIC = 16 μg/mL is intermediate and a MIC ≥ 32 μg/mL is resistant. ⁹Also consistent with other *Enterobacter* species

Figure 1: Colony Morphology



Date: 01 SEP 2016

Signature:

BEI Resources Authentication

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org Tel: 800-359-7370

¹⁰Purity of this lot was assessed for 7 days at 37°C in an aerobic atmosphere on Tryptic Soy agar.