SUPPORTING INFECTIOUS DISEASE RESEARCH

### Mycobacterium tuberculosis, Strain 11895-0

#### Catalog No. NR-48755

**Product Description:** *Mycobacterium tuberculosis* (*M. tuberculosis*), strain 11895-0 was isolated in October 2012 from a subculture of a strain originally isolated from a patient with pulmonary tuberculosis in the Republic of South Africa. *M. tuberculosis*, strain 11895-0 was deposited as a multidrug-resistant (MDR) strain with resistance to amikacin, capreomycin, ethambutol, ethionamide, isoniazid, kanamycin, pyrazinamide, rifampin and streptomycin.

#### Lot<sup>1</sup>: 63950984

## Manufacturing Date: 15JAN2016

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis <sup>2</sup>		
Cellular morphology	Gram-positive rods	Gram-positive rods
Colony morphology <sup>3</sup>	Report results	Irregular, raised, entire, rough and
e clony melphology		cream
Growth rate	≥ 7 days	22 days
Growth at 26°C	Negative	Negative
Growth at 37°C	Positive	Positive
Acid-fast stain	Positive (red colonies)	Positive (red colonies)
Pigmentation in the dark (Scotochromogen)	Negative (no pigment)	Negative (no pigment)
Photoinduction for 1 hour (Photochromogen)	Negative (no pigment)	Negative (no pigment)
Nonchromogen (no pigment)	Positive (no pigment)	Positive (no pigment)
Biochemical tests		
Niacin production <sup>4</sup>	Positive	Positive
Nitrate reduction	Positive	Positive
Pyrazinamidase	Positive	Positive
•		
Antibiotic Susceptibility Profile		
Sensititre <sup>™</sup> System <sup>5,6</sup>		
Amikacin	Report results	> 16 µg/mL
Cycloserine	Report results	32 µg/mL
Ethambutol	Report results	8 μg/mL <sup>7</sup>
Ethionamide	Report results	5 μg/mL <sup>7</sup>
Isoniazid	Report results	> 4 µg/mL
Kanamycin	Report results	> 40 µg/mL
Moxifloxacin	Report results	0.5 μg/mL
Ofloxacin	Report results	4 μg/mL
Para-aminosalicylic acid	Report results	1 μg/mL <sup>7</sup>
Rifabutin	Report results	0.5 μg/mL <sup>7</sup>
Rifampin	Report results	> 16 µg/mL
Streptomycin	Report results	> 32 µg/mL <sup>7</sup>
Constunia Analysia		
Genotypic Analysis		100% appropriate identity to
Sequencing of Heat Shock Protein 65 gene	$\geq$ 99% sequence identity to	100% sequence identity to
(~ 310 base pairs)	<i>M. tuberculosis</i> type strain	<i>M. tuberculosis</i> type strain
	(GenBank: AL123456)	(GenBank: AL123456) <sup>8</sup>
Purity (post-freeze)		
Middlebrook 7H10 agar with OADC enrichment <sup>9</sup>	Growth consistent with expected	Growth consistent with expected
madiobrook intro agar with onbo childhillent	colony morphology	colony morphology
Tryptic Soy agar <sup>9</sup>	Report results	No growth
Viability (post-freeze) <sup>3</sup>	Growth	Growth

<sup>1</sup>NR-48755 was produced by inoculation of the deposited material into Middlebrook 7H9 broth with ADC enrichment. Broth inoculum was added to Middlebrook 7H10 agar with OADC enrichment kolles, which were grown for 29 days at 37°C in an aerobic atmosphere with 5% CO<sub>2</sub> to produce this lot.

BEI Resources www.beiresources.org E-mail: <u>contact@beiresources.org</u> Tel: 800-359-7370 Fax: 703-365-2898 dieii resources

# **Certificate of Analysis for NR-48755**

SUPPORTING INFECTIOUS DISEASE RESEARCH

<sup>2</sup>Information on Mycobacterium testing is available from Ribón, W. "Biochemical Isolation and Identification of Mycobacteria." <u>Biochemical Testing</u>. (2012) Jose C. Jimenez-Lopez (Ed.), InTech, <u>http://www.intechopen.com/books/biochemical-testing/biochemical-isolation-and-identification-of-mycobacteria</u> and Lévy-Frébault, V. V. and F. Portaels. "Proposed Minimal Standards for the Genus *Mycobacterium* and for Description of New Slowly Growing *Mycobacterium* Species." <u>Int. J. Syst. Bacteriol.</u> 42 (1992): 315-323. PubMed: 1581193.

<sup>3</sup>22 days at 37°C in an aerobic atmosphere with 5% CO<sub>2</sub> on Middlebrook 7H10 agar with OADC enrichment

<sup>4</sup>All mycobacteria produce niacin but only *M. tuberculosis* accumulates it, resulting in a positive test for *M. tuberculosis*.

<sup>5</sup>Sensititre<sup>TM</sup> System *Mycobacterium tuberculosis* MIC Plate, Thermo Scientific<sup>TM</sup>, catalog number MYCOTB

<sup>6</sup>Minimum Inhibitory Concentration (MIC); No Clinical & Laboratory Standards Institute (CLSI) interpretations of the Sensititre<sup>™</sup> System data for *M. tuberculosis* are currently available.

<sup>7</sup>For streptomycin, ethionamide, para-aminosalicylic acid, rifabutin and ethambutol, the endpoint for these drugs is determined by the well with approximately 80% inhibition of growth compared to the positive control well with no drug.

<sup>8</sup>Also consistent with *M. africanum*, *M. bovis*, *M. canettii*, *M. caprae* and *M. microti* 

<sup>9</sup>Purity of this lot was assessed for 22 days at 37°C in an aerobic atmosphere with 5% CO<sub>2</sub>.

Date: 10 AUG 2017

Signature:

**BEI Resources Authentication** 

ATCC<sup>®</sup>, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC<sup>®</sup>'s knowledge.

ATCC<sup>®</sup> is a trademark of the American Type Culture Collection. You are authorized to use this product for research use only. It is not intended for human use.

