

# **Product Information Sheet for NR-4617**

Genomic DNA from Salmonella enterica subsp. enterica, 2004 Pennsylvania Tomato Outbreak, Serovar Typhimurium, Isolate 4

Catalog No. NR-4617

# For research use only. Not for human use.

#### **Contributor:**

Carol H. Sandt, Molecular Microbiology Section, Bureau of Laboratories, Pennsylvania Department of Health, Lionville, Pennsylvania

### **Product Description:**

Genomic DNA was obtained from a preparation of *Salmonella enterica* (*S. enterica*) subsp. *enterica* serovar Typhimurium that was isolated from the stool of a patient with diarrhea during the 2004 *Salmonella* outbreak in Pennsylvania.<sup>1,2</sup>

S. enterica subsp. enterica serovar Typhimurium is a major cause of gastroenteritis. These bacteria are host generalists that occur in humans and many other mammals. Additionally, this serovar causes typhoid-like disease in mice and is used as a mouse model of human typhoid fever.<sup>3</sup> The complete genome sequence of several strains of S. enterica subsp. enterica serovar Typhimurium is in progress [strain DT104 (Definitive Type 104; a multidrug resistant strain), strain SL1344 (a genetically marked subline of a calf-virulent isolate) and strain TR7095 (a wild-type strain)] and strain LT2 has been completed (GenBank: AE006468).<sup>3</sup>

NR-4617 has been qualified for PCR applications by amplification of approximately 1500 bp of the 16S ribosomal RNA.

#### **Material Provided:**

Each vial contains 4 to 6  $\mu$ g of bacterial genomic DNA in TE buffer (10 mM Tris-HCl and 1 mM EDTA, pH 7.4). The concentration is shown on the Certificate of Analysis. The vial should be centrifuged prior to opening.

## Packaging/Storage:

NR-4617 was packaged aseptically in screw-capped plastic cryovials. The product is provided frozen on dry ice and should be stored at -20°C or colder immediately upon arrival. Freeze-thaw cycles should be minimized.

### Citation:

Acknowledgment for publications should read "The following reagent was obtained through the NIH Biodefense and Emerging Infections Research Resources Repository, NIAID, NIH: Genomic DNA from Salmonella enterica subsp. enterica, 2004 Pennsylvania Tomato Outbreak, Serovar

Typhimurium, Isolate 4, NR-4617."

## Biosafety Level: 1

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the following publication: U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. Biosafety in Microbiological and Biomedical Laboratories. 5th ed. Washington, DC: U.S. Government Printing Office, 2007; see www.cdc.gov/od/ohs/biosfty/bmbl5/bmbl5toc.htm.

#### **Disclaimers:**

You are authorized to use this product for research use only. It is not intended for human use.

Use of this product is subject to the terms and conditions of the BEI Resources Material Transfer Agreement (MTA). The MTA is available on our Web site at <a href="https://www.beiresources.org">www.beiresources.org</a>.

While BEI Resources uses reasonable efforts to include accurate and up-to-date information on this product sheet, neither ATCC® nor the U.S. Government make any warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. Neither ATCC® nor the U.S. Government warrants that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, use and disposal. ATCC® and the U.S. Government are not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, the U.S. Government, ATCC®, their suppliers and contributors to BEI Resources are not liable for damages arising from the misidentification or misrepresentation of products.

#### **Use Restrictions:**

This material is distributed for internal research, non-commercial purposes only. This material, its product or its derivatives may not be distributed to third parties. Except as performed under a U.S. Government contract, individuals contemplating commercial use of the material, its products or its derivatives must contact the contributor to determine if a license is required. U.S. Government contractors may need a license before first commercial sale.

### References:

- Sandt, C. H., et al. "The Key Role of Pulsed-Field Gel Electrophoresis in Investigation of a Large Multiserotype and Multistate Food-Borne Outbreak of Salmonella Infections Centered in Pennsylvania." J. Clin. Microbiol. 44 (2006): 3208-3212. PubMed: 16954249.
- Centers for Disease Control and Prevention (CDC). "Outbreaks of Salmonella Infections Associated with

Biodefense and Emerging Infections Research Resources Repository

P.O. Box 4137 Manassas, VA 20108-4137 USA www.beiresources.org Fax: 703-365-2898 E-mail: contact@beiresources.org

© 2008 American Type Culture Collection (ATCC). All rights reserved. Page 1 of 2

800-359-7370



# **Product Information Sheet for NR-4617**

- Eating Roma Tomatoes--United States and Canada, 2004." Morb. Mortal. Wkly. Rep. 54 (2005): 325-328. PubMed: 15815562.
- McClelland, M., et al. "Complete Genome Sequence of Salmonella enterica Serovar Typhimurium LT2." <u>Nature</u> 413 (2001): 852–856. PubMed: 11677609. GenBank: AE006468.

ATCC® is a trademark of the American Type Culture Collection.

800-359-7370

Fax: 703-365-2898

E-mail: contact@beiresources.org