

Kilbourne F7: A/New Jersey/11/76 (H1N1) Mutant, High (H) Yield

Catalog No. NR-3598

Derived from NIAID Catalog No. V-331-0E5462

For research use only. Not for human use.

Contributor:

National Institutes of Allergy and Infectious Diseases,
National Institutes of Health

Product Description:

Virus Classification: *Orthomyxoviridae, Influenzavirus A*

Species: Influenza A virus

Mutant: A/New Jersey/11/76 (H1N1) Mutant, High (H) Yield (Kilbourne F7)¹⁻³

Comments: This is a high (H) yield HA mutant of A/New Jersey/11/76 (H1N1) virus that was isolated from the 1976 epidemic at Ft. Dix, New Jersey.⁴ This strain also grows to high (hy) yield in chick embryos.

Material Provided:

Each vial contains approximately 1 mL of pooled allantoic fluid from specific-pathogen free (SPF) embryonated chicken eggs infected with recombinant influenza A virus mutant of A/New Jersey/11/76 (H1N1).

Packaging/Storage:

NR-3598 was packaged aseptically in screw-capped plastic cryovials. The product is provided frozen and should be stored at -70°C or colder immediately upon arrival. For long-term storage, the vapor phase of a liquid nitrogen freezer is recommended. Freeze-thaw cycles should be avoided.

Growth Conditions:

Host: 9 to 11-day-old SPF embryonated chicken eggs

Infection: Embryonated chicken eggs must be candled for viability prior to inoculation

Incubation: 1-3 days at 33-35°C in a humidified chamber without CO₂

Effect: Hemagglutination activity using chicken red blood cells and allantoic fluid from infected embryonated chicken eggs

Citation:

Acknowledgment for publications should read "The following reagent was obtained through the NIH Biodefense and Emerging Infections Research Resources Repository, NIAID, NIH: Kilbourne F7: A/New Jersey/11/76 (H1N1) Mutant, High (H) Yield, NR-3598."

Biosafety Level: 2

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the following publication: U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and

Prevention, and National Institutes of Health. Biosafety in Microbiological and Biomedical Laboratories. 5th ed. Washington, DC: U.S. Government Printing Office, 2007; see www.cdc.gov/od/ohs/biosfty/bmb15/bmb15toc.htm.

Disclaimers:

You are authorized to use this product for research use only. It is not intended for human use.

Use of this product is subject to the terms and conditions of the BEI Resources Material Transfer Agreement (MTA). The MTA is available on our Web site at www.beiresources.org.

While BEI Resources uses reasonable efforts to include accurate and up-to-date information on this product sheet, neither ATCC® nor the U.S. Government make any warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. Neither ATCC® nor the U.S. Government warrants that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, use and disposal. ATCC® and the U.S. Government are not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, the U.S. Government, ATCC®, their suppliers and contributors to BEI Resources are not liable for damages arising from the misidentification or misrepresentation of products.

Use Restrictions:

This material is distributed for internal research, non-commercial purposes only. This material, its product or its derivatives may not be distributed to third parties. Except as performed under a U.S. Government contract, individuals contemplating commercial use of the material, its products or its derivatives must contact the contributor to determine if a license is required. U.S. Government contractors may need a license before first commercial sale.

References:

1. http://www.flu-archive.org/data_sheets/F7.doc
2. <http://www.flu-archive.org/>
3. http://www.flu-archive.org/search/results.pl?search_string=&join_type=and
4. Kilbourne, E. D. "Genetic Dimorphism in Influenza Viruses: Characterization of Stably Associated Hemagglutinin Mutants Differing in Antigenicity and Biological Properties." *Proc. Natl. Acad. Sci. U. S. A.* 75 (1978): 6258-6262. PubMed: 282644.
5. Gambaryan, A. S., et al. "Differences in the Biological Phenotype of Low-Yielding (L) and High-Yielding (H) Variants of Swine Influenza Virus A/NJ/11/76 are Associated with Their Different Receptor-Binding Activity." *Virology* 247 (1998): 223-231. PubMed:

- 9705915.
6. Kilbourne E. D., B. C. Easterday and S. McGregor. "Evolution to Predominance of Swine Influenza Virus Hemagglutinin Mutants of Predictable Phenotype during Single Infections of the Natural Host." Proc. Natl. Acad. Sci. U. S. A. 85 (1988): 8098-8101. PubMed: 3186713.
 7. Kilbourne E. D., et al. "Hemagglutinin Polymorphism as the Basis for Low- and High-Yield Phenotypes of Swine Influenza Virus." Proc. Natl. Acad. Sci. U. S. A. 85 (1988): 7782-7785. PubMed: 3174662.
 8. Kilbourne E. D., G. W. Both and W. Gerhard. "Pleiotropic Effects of a Single Amino Acid Change on Antigenicity and Biologic Function of Swine Influenza Virus Hemagglutinin Mutants." In: Negative Strand Viruses. Eds. D. H. Bishop and R. W. Compans. Academic Press New York, 1984. 233-237.
 9. Kilbourne, E. D. "The Genetic Dissection of Viral Virulence." Ann N. Y. Acad. Sci. 435 (1984): 32-38.
 10. Kilbourne E. D., W. Gerhard and C. W. Whitaker. "Monoclonal Antibodies to the Hemagglutinin Sa Antigenic Site of A/PR/8/34 Influenza Virus Distinguish Biologic Mutants of Swine Influenza Virus." Proc. Natl. Acad. Sci. U. S. A. 80 (1983): 6399-6402. PubMed: 6194531.
 11. Kilbourne E. D., S. McGregor and B. C. Easterday. "Transmission in Swine of Hemagglutinin Mutants of Swine Influenza Virus." In: The Replication of Negative Strand Viruses. Eds. D. H. Bishop and R. W. Compans. Elsevier North Holland, Inc., 1981. 449-453.
 12. Kilbourne, E. D. "Influenza: Viral Determinants of the Pathogenicity and Epidemicity of an Invariant Disease of Variable Occurrence." Philos. Trans. R. Soc. Lond. B. Biol. Sci. 288 (1980): 291-297. PubMed: 6103545.
 13. Kilbourne E. D., S. McGregor and B. C. Easterday. "Hemagglutinin Mutants of Swine Influenza Virus Differing in Replication Characteristics in their Natural Host." Infect. Immun. 26 (1979): 197-201. PubMed: 227792.
 14. Kilbourne, E. D. "Molecular Epidemiology-Influenza as Archetype." Harvey Lect. 73 (1979): 225-258. PubMed: 396276.

ATCC® is a trademark of the American Type Culture Collection.

