

Certificate of Analysis for NR-43756

Genomic RNA from Influenza A Virus, A/Uruguay/716/2007 (H3N2)

Catalog No. NR-43756


Product Description: Genomic RNA was extracted from a preparation of pooled allantoic fluid from specific pathogen free embryonated chicken eggs¹ infected with influenza A virus, A/Uruguay/716/2007 (H3N2).

Lot²: 70015822 Manufacturing Date: 01JUN2018

TEST	SPECIFICATIONS	RESULTS
Sequencing of Influenza A Specific Region		
Matrix gene (~ 930 base pairs)	≥ 98% identity with A/Uruguay/716/2007 (H3N2) (GenBank: EU716434) Consistent with source organism (BEI Resources NR-42003)	100% identity with A/Uruguay/716/2007 (H3N2) (GenBank: EU716434) Consistent with NR-42003
Functional Activity by RT-PCR Amplification ³		
Matrix gene	~ 1030 base pair amplicon	~ 1030 base pair amplicon (Figure 1)
Pre-Vial Concentration by RiboGreen® Measurement (Viral, Cellular and Carrier)4	Report results	99 ng per 100 μL (0.99 μg/mL)
Estimated Amount per Vial ⁴	Report Results	99 ng
Virus Inactivation 10% of total yield incubated in SPF embryonated chicken eggs ¹ and evaluated by hemagglutination assay ⁵	No viable virus detected	No viable virus detected
Mycoplasma Contamination DNA detection by PCR of extracted Test Article	None detected	None detected

¹10-day-old SPF Embryonated Chicken Eggs acquired from Charles River, Norwich, Connecticut, USA

Figure 1: Functional Activity of NR-43756 by RT-PCR Amplification of Matrix gene



Lane 1: Invitrogen™ TrackIt™ 1 Kb Plus DNA Ladder

Lane 2: PCR product from 1 µL of NR-43756

BEI Resources www.beiresources.org E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

²Nucleic acid was extracted from a preparation of influenza A virus, Uruguay/716/2007 (H3N2) (BEI Resources NR-42003 lot 70010871) using a QIAamp[®] Viral RNA Mini Kit (Qiagen 52906).

³Amplified using iTaq™ Universal SYBR® Green One-Step Kit (Bio-Rad 172-5151) with 5 µL of NR-43756 in a 50 µL reaction

⁴Measurement is determined pre-vial prior to dilution due to the limit of detection of the quantification method.

⁵Use of the QIAamp[®] Viral RNA Mini Kit has been demonstrated to consistently inactivate 100% of influenza viruses.



Certificate of Analysis for NR-43756

/Heather Couch/

Heather Couch 10 SEP 2018

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

Provided by NIAID

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898

NR-43756 70015822 10SEP2018