

**Genomic DNA from *Mycobacterium africanum*, Strain NLA009502090**

**Catalog No. NR-49655**

**Product Description:**

Genomic DNA was extracted from a preparation of *Mycobacterium africanum* (*M. africanum*), strain NLA009502090.

**Lot: 70003664<sup>1,2</sup>**

**Manufacturing Date: 07NOV2017**

| TEST   | SPECIFICATIONS   | RESULTS   |
|--|--|---|
| <b>Genotypic Analysis</b><br>Sequencing of 16S ribosomal RNA gene<br>(~ 1410 base pairs) | ≥ 99% sequence identity to<br><i>M. africanum</i> type strain<br>(GenBank: AF480605.1) | 99.9% sequence identity to<br><i>M. africanum</i> type strain<br>(GenBank: AF480605.1) <sup>3</sup> |
| <b>Agarose Gel Electrophoresis</b>   | High molecular weight chromosomal<br>DNA   | High molecular weight chromosomal<br>DNA (Figure 1)   |
| <b>Concentration by PicoGreen® Measurement</b>   | 0.7 to 1.5 µg in 25 to 100 µL per vial   | 0.2 µg in 38 µL per vial (6 µg/mL)  |
| <b>Amount per Vial</b>   | 0.7 to 1.5 µg  | 0.2 µg <sup>4</sup>   |
| <b>Functional Activity by PCR Amplification</b><br>16S ribosomal RNA gene                | ~ 1500 base pair amplicon  | ~ 1500 base pair amplicon   |
| <b>OD<sub>260</sub>/OD<sub>280</sub> Ratio</b>   | 1.7 to 2.1   | 1.5 <sup>5</sup>  |
| <b>Bacterial Inactivation</b><br>10% of total yield plated on agar <sup>6,7</sup>        | No viable bacteria detected  | No viable bacteria detected   |

<sup>1</sup>The bacterial preparation used for extraction of genomic DNA was produced from the BEI Resources NRS-49261 lot 70003658. Genomic DNA was extracted using proprietary technology.

<sup>2</sup>NR-49655 lot 70003664 was vialled in TE buffer (10 mM Tris-HCl, 1 mM EDTA, pH ~ 8.0).

<sup>3</sup>Also consistent with other *Mycobacterium* species


<sup>4</sup>The amount of genomic DNA in the vial falls below the current specifications, but does not negatively impact the final product.

<sup>5</sup>Although the OD<sub>260</sub>/OD<sub>280</sub> ratio falls below the current specification, the material was still functional in PCR applications.

<sup>6</sup>37 days at 37°C in an aerobic atmosphere with 5% CO<sub>2</sub> on Middlebrook 7H10 agar with OADC enrichment.

<sup>7</sup>An extraction procedure was used that has been shown to consistently inactivate 100% of Gram-positive and Gram-negative bacteria.

**Figure 1: Agarose Gel Electrophoresis**



Lane 1: Invitrogen™ TrackIt™ 100 bp DNA Ladder  
Lane 2: 120 ng of NR-49655

/Heather Couch/

Heather Couch

04 SEP 2019

Program Manager or designee, ATCC Federal Solutions

ATCC<sup>®</sup>, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC<sup>®</sup>'s knowledge.

ATCC<sup>®</sup> is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

