

Certificate of Analysis for NR-51611

Pseudomonas aeruginosa, Strain MRSN 409937

Catalog No. NR-51611

This reagent is the tangible property of the U.S. Government.

Product Description:

Pseudomonas aeruginosa (P. aeruginosa), strain MRSN 409937 was isolated in 2016 from human fluid as part of a surveillance program in the United States. P. aeruginosa, strain MRSN 409937 was deposited as sensitive to amikacin, gentamicin, imipenem and tobramycin and resistant to aztreonam, cefepime, ceftazidime, ciprofloxacin, levofloxacin, meropenem and piperacillin/tazobactam.

Lot: 70025128¹ Manufacturing Date: 26JUL2019

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis		
Cellular morphology	Gram-negative rods	Gram-negative rods
Colony morphology ²	Report results	Circular, flat, undulate, opaque and green (Figure 1)
Motility (wet mount)	Report results	Motile
VITEK® 2 (GN card)	P. aeruginosa (≥ 89%)	P. aeruginosa (98%)
Antibiotic Susceptibility Profile ³ VITEK® (AST-GN81 Card)		
Ampicillin	Report results	Resistant (≥ 32 µg/mL)
Amoxicillin/clavulanic Acid	Report results	Resistant (≥ 32 µg/mL)
Piperacillin/tazobactam	Resistant	Resistant (≥ 128 µg/mL)
Cefazolin	Report results	Resistant (≥ 64 µg/mL)
Cefoxitin	Report results	Resistant (≥ 64 µg/mL)
Ceftazidime	Resistant	Resistant (≥ 64 µg/mL)
Ceftriaxone	Report results	Resistant (≥ 64 µg/mL)
Cefepime	Resistant	Resistant (≥ 64 µg/mL)
Meropenem	Resistant	Intermediate (4 µg/mL) ⁴
Amikacin	Sensitive	Sensitive (4 µg/mL)
Gentamicin	Sensitive	Sensitive (≤ 1 µg/mL)
Tobramycin	Sensitive	Sensitive (≤ 1 µg/mL)
Ciprofloxacin	Resistant	Resistant (≥ 4 µg/mL)
Levofloxacin	Resistant	Resistant (≥ 8 µg/mL)
Tetracycline	Report results	Resistant (≥ 16 µg/mL)
Nitrofurantoin	Report results	Resistant (≥ 512 µg/mL)
Trimethoprim/sulfamethoxazole	Report results	≥ 320 µg/mL ⁵
Genotypic Analysis Sequencing of 16S ribosomal RNA gene (~ 1420 base pairs)	≥ 99% sequence identity to P. aeruginosa, strain MRSN 409937 (GenBank: RXTX01000079.1)	99.9% sequence identity to P. aeruginosa, strain MRSN 409937 (GenBank: RXTX01000079.1)
Purity (post-freeze) ⁶	Growth consistent with expected colony morphology	Growth consistent with expected colony morphology
Viability (post-freeze) ²	Growth	Growth
	<u> </u>	

¹NR-51611 was produced by inoculation of the depositor material into Tryptic Soy broth and grown for 1 day at 37°C in an aerobic atmosphere. Broth inoculum was added to Tryptic Soy agar kolles, which were grown for 1 day at 37°C in an aerobic atmosphere to produce this lot.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org Tel: 800-359-7370

Fax: 703-365-2898

²1 day at 37°C in an aerobic atmosphere on Tryptic Soy agar

³Minimum Inhibitory Concentration (MIC); MIC Interpretation Guideline: CLSI M100-S28 (2018)

⁴P. aeruginosa, strain MRSN 409937 was deposited as resistant to meropenem. Antibiotic susceptibility testing performed in duplicate identified strain MRSN 409937 as having an intermediate resistance to meropenem.

⁵Trimethoprim/sulfamethoxazole MIC interpretive standards are not available for *P. aeruginosa*, however most clinical isolates are resistant to trimethoprim/sulfamethoxazole. For more information, please refer to Köhler, T., et al. "Multidrug Efflux in Intrinsic Resistance to Trimethoprim and Sulfamethoxazole in *Pseudomonas aeruginosa*." <u>Antimicrob. Agents Chemother.</u> 40 (1996): 2288-2290. PubMed: 9036831.



Certificate of Analysis for NR-51611

⁶Purity of this lot was assessed for 7 days at 37°C in an aerobic atmosphere with and without 5% CO₂ on Tryptic Soy agar.

Figure 1: Colony Morphology



/Heather Couch/ Heather Couch

12 DEC 2019

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

BEI Resources www.beiresources.org E-mail: contact@beiresources.org
Tel: 800-359-7370

Fax: 703-365-2898