

***Candida albicans*, Strain P78048**

Catalog No. NR-29434

Product Description:

Candida albicans (*C. albicans*), strain P78048 was isolated in 2000 from the bloodstream of a human in Winnipeg, Manitoba, Canada. *C. albicans*, strain P78048 is reported to be a member of genetic clade I with an alpha/alpha MTL genotype. NR-29434 was produced by inoculation of BEI Resources NR-29434 lot 60842872 onto Yeast Mold agar, which was grown for 2 days at 25°C in an aerobic atmosphere. The agar growth was harvested with 10% glycerol to produce this lot.

Lot: 70016187

Manufacturing Date: 05JUL2018

TEST	SPECIFICATIONS	RESULTS
Phenotypic Analysis Cellular morphology (pre-freeze) 2 days at 25°C in an aerobic atmosphere on Yeast Mold agar Colony morphology (pre-freeze) CHROMagar™ Candida Testing ¹ Biochemical Tests VITEK® 2 (YST card) VITEK® MS (MALDI-TOF)	Report results Report results Green colonies (<i>C. albicans</i>) <i>C. albicans</i> (≥ 89%) <i>C. albicans</i>	Subglobose to ovoid, in singles and pairs; no pseudohyphae observed (Figure 1) Butyrous, dull, matte and off-white (Figure 2) Green colonies (<i>C. albicans</i>) <i>C. albicans</i> (99%) <i>C. albicans</i> (99.9%)
Genotypic Analysis Sequencing of partial 18S ribosomal RNA (rRNA) gene, internal transcribed spacer (ITS) 1, 5.8S rRNA gene, ITS 2, partial 26S rRNA (~ 510 base pairs) Sequencing of 28S rRNA gene (~ 610 base pairs)	≥ 99% sequence identity to <i>C. albicans</i> , strain SC5314 (GenBank: CP017630.1) ≥ 99% sequence identity to <i>C. albicans</i> , strain SC5314 (GenBank: CP017630.1)	100% sequence identity to <i>C. albicans</i> , strain SC5314 (GenBank: CP017630.1) 100% sequence identity to <i>C. albicans</i> , strain SC5314 (GenBank: CP017630.1) ²
Purity³ Nutrient broth with 0.1% Yeast Extract at 25°C Nutrient broth with 0.1% Yeast Extract at 37°C	No bacterial growth No bacterial growth	Non-turbid growth Non-turbid growth
Viability (post-freeze)	Growth	Growth

¹CHROMagar™ testing was performed on BEI Resources NR-29434 lot 60842872. CHROMagar™ displays high sensitivity and specificity for three of the major *Candida* species (Odds, F. C., and R. Bernarts. "CHROMagar Candida, a New Differential Isolation Medium for the Presumptive Identification of Clinically Important *Candida* species." *J. Clin. Microbiol.* 32 (1994): 1923-1929. PubMed: 7989544.).

²Also consistent with *C. africana*

³Clarity of broth was determined by visual inspection after 2 days in an aerobic atmosphere.

Figure 1: Cellular Morphology




Figure 2: Colony Morphology



/Heather Couch/
Heather Couch

10 MAR 2020

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

