

***Borrelia miyamotoi*, Strain HT31**

**Catalog No. NR-51675**

**Product Description:**

*Borrelia miyamotoi* (*B. miyamotoi*), strain HT31 was isolated from the abdomen of an unfed female *Ixodes persulcatus* tick collected from vegetation between 1990 and 1992 in Shiretoko, Hokkaido, Japan. The deposited material was inoculated into Revised Barbour-Stoenner-Kelly broth and grown for two passages at 33°C in an aerobic atmosphere with 5% CO<sub>2</sub>, and the resulting subculture was vialled and frozen. NR-51675 lot 70032899 was produced by inoculation of the deposited material into Revised Barbour-Stoenner-Kelly broth. After two passages, the culture was grown for 7 days at 33°C in an aerobic atmosphere with 5% CO<sub>2</sub>. Broth inoculum was added to Revised Barbour-Stoenner-Kelly broth and grown for 6 days at 33°C in a microaerophilic atmosphere to produce this lot.


**Lot: 70032899**

**Manufacturing Date: 29JAN2020**

TEST	SPECIFICATIONS	RESULTS
<b>Phenotypic Analysis</b> Cellular morphology 7 days at 33°C in a microaerophilic atmosphere in Revised Barbour-Stoenner-Kelly broth Motility (wet mount)	Spirochete  Report results	Spirochete  Motile
<b>Genotypic Analysis</b> Sequencing of 16S ribosomal RNA (rRNA) gene (1410 base pairs)	≥ 99% sequence identity to <i>B. miyamotoi</i> , strain HT31 (GenBank: AB904793.1)	100% sequence identity to <i>B. miyamotoi</i> , strain HT31 (GenBank: AB904793.1)
<b>Purity (post-freeze)</b> 7 days at 33°C in an aerobic atmosphere with 5% CO <sub>2</sub> on Tryptic Soy agar with 5% defibrinated sheep blood 7 days at 37°C in an aerobic atmosphere with 5% CO <sub>2</sub> on Tryptic Soy agar with 5% defibrinated sheep blood	Growth consistent with colony morphology or no growth No growth	No growth No growth
<b>Viability (post-freeze)</b> Visual observation 7 days at 33°C in a microaerophilic atmosphere in Revised Barbour-Stoenner-Kelly broth LIVE/DEAD® BacLight™ Bacterial Viability	Growth  Green fluorescence visible	Growth  Green fluorescence visible (Figure 1) <sup>1</sup>

<sup>1</sup>Determined after 7 days at 33°C in an aerobic atmosphere with 5% CO<sub>2</sub> in Revised Barbour-Stoenner-Kelly broth with LIVE/DEAD® BacLight™ Bacterial Viability Kit, 1000× magnification (Invitrogen™ L34856). Cells with a compromised membrane that are dead or dying will stain red, while cells with an intact membrane will stain green.

Figure 1: LIVE/DEAD® BacLight™ Bacterial Viability



/Heather Couch/  
Heather Couch

28 APR 2020

Program Manager or designee, ATCC Federal Solutions

ATCC®, on behalf of BEI Resources, hereby represents and warrants that the material provided under this certificate has been subjected to the tests and procedures specified and that the results described, along with any other data provided in this certificate, are true and accurate to the best of ATCC®'s knowledge.

ATCC® is a trademark of the American Type Culture Collection.

You are authorized to use this product for research use only. It is not intended for human use.

