

# **Product Information Sheet for NR-44267**

## Mycobacterium intracellulare, Strain 1956

## Catalog No. NR-44267

### For research use only. Not for human use.

### Contributor:

Diane Ordway, Ph.D., Assistant Professor, Department of Microbiology, Immunology and Pathology, Colorado State University, Fort Collins, Colorado, USA and National Institute of Allergy and Infectious Diseases (NIAID), National Institutes of Health (NIH), Bethesda, Maryland, USA

#### Manufacturer:

**BEI Resources** 

### **Product Description:**

Bacteria Classification: Mycobacteriaceae, Mycobacterium

Species: Mycobacterium intracellulare

<u>Strain</u>: 1956

Original Source: Mycobacterium intracellulare (M. intracellulare), strain 1956 was isolated in 2011 from human sputum at NIAID, NIH, Bethesda, Maryland, USA.<sup>1</sup>

<u>Comment</u>: *M. intracellulare*, strain 1956 is part of the Top Priority Nontuberculous Mycobacteria Whole Genome Sequencing Project at the <u>Genomic Sequencing Center for Infectious Diseases</u> (GSCID) at University of Maryland School of Medicine. The complete genome sequence of *M. intracellulare*, strain 1956 has been sequenced (GenBank: JAOG00000000).

*M. intracellulare* is an acid-fast, Gram-positive, non-motile, non-pigmenting, slow-growing bacillus ubiquitous in the environment and typically found in water, soil and hospital wards.<sup>2</sup> It is frequently associated with human pulmonary infections in immunocompromised (often HIV) patients. Along with *Mycobacterium avium* (*M. avium*), *M. intracellulare* is a member of the *M. avium* complex (MAC), the nontuberculous mycobacteria most often isolated in clinical settings.<sup>2</sup> Differentiating between *M. intracellulare* and *M. avium* is often challenging, as few distinct disparities exist between them.<sup>3,4,5,6</sup> Given the difficulty and time to identify these two similar species, clinical reports commonly cite them together as *M. avium-intracellulare* (MAI) or MAC infections.<sup>7</sup>

### **Material Provided:**

Each vial contains approximately 0.5 mL of bacterial culture in Middlebrook 7H9 broth with ADC Enrichment supplemented with 10% glycerol.

<u>Note</u>: If homogeneity is required for your intended use, please purify prior to initiating work.

### Packaging/Storage:

NR-44267 was packaged aseptically in screw-capped plastic cryovials. The product is provided frozen and should be stored at -60°C or colder immediately upon arrival. For long-term storage, the vapor phase of a liquid nitrogen freezer is recommended. Freeze-thaw cycles should be avoided.

#### **Growth Conditions:**

Media

Middlebrook 7H9 broth with ADC enrichment or equivalent Middlebrook 7H10 agar with OADC enrichment or Lowenstein-Jensen agar or equivalent

Incubation:

Temperature: 37°C

Atmosphere: Aerobic with 5% CO<sub>2</sub>

Propagation:

- 1. Keep vial frozen until ready for use; then thaw.
- Transfer the entire thawed aliquot into a single tube of broth.
- Use several drops of the suspension to inoculate an agar slant and/or plate.
- Incubate the tube, slant and/or plate at 37°C for 2 to 8 weeks.

#### Citation:

Acknowledgment for publications should read "The following reagent was obtained through BEI Resources, NIAID, NIH: *Mycobacterium intracellulare*, Strain 1956, NR-44267."

### Biosafety Level: 2

Appropriate safety procedures should always be used with this material. Laboratory safety is discussed in the following publication: U.S. Department of Health and Human Services, Public Health Service, Centers for Disease Control and Prevention, and National Institutes of Health. Biosafety in Microbiological and Biomedical Laboratories. 5th ed. Washington, DC: U.S. Government Printing Office, 2009; see www.cdc.gov/biosafety/publications/bmbl5/index.htm.

This publication recommends that practices with this agent include the use of respiratory protection and the implementation of specific procedures and use of specialized equipment to prevent and contain aerosols.

### **Disclaimers:**

You are authorized to use this product for research use only. It is not intended for human use.

Use of this product is subject to the terms and conditions of the BEI Resources Material Transfer Agreement (MTA). The MTA is available on our Web site at www.beiresources.org.

While BEI Resources uses reasonable efforts to include accurate and up-to-date information on this product sheet, neither ATCC® nor the U.S. Government makes any warranties or representations as to its accuracy. Citations from scientific literature and patents are provided for informational purposes only. Neither ATCC® nor the U.S. Government warrants that such information has been confirmed to be accurate.

This product is sent with the condition that you are responsible for its safe storage, handling, use and disposal. ATCC® and the U.S. Government are not liable for any damages or injuries arising from receipt and/or use of this product. While reasonable effort is made to ensure authenticity and reliability of materials on deposit, the U.S. Government, ATCC®, their

BEI Resources www.beiresources.org E-mail: contact@beiresources.org

Tel: 800-359-7370 Fax: 703-365-2898

NR-44267 22MAY2020



# **Product Information Sheet for NR-44267**

suppliers and contributors to BEI Resources are not liable for damages arising from the misidentification or misrepresentation of products.

### **Use Restrictions:**

This material is distributed for internal research, non-commercial purposes only. This material, its product or its derivatives may not be distributed to third parties. Except as performed under a U.S. Government contract, individuals contemplating commercial use of the material, its products or its derivatives must contact the contributor to determine if a license is required. U.S. Government contractors may need a license before first commercial sale.

### References:

- 1. Ordway, D., Personal Communication.
- Inderlied, C. B., C. A. Kemper and L. E. Bermudez. "The Mycobacterium avium Complex." <u>Clin. Microbiol. Rev.</u> 6 (1993): 266-310. PubMed: 8358707.
- Renvoisé, A., et al. "Significant Difference in Drug Susceptibility Distribution between Mycobacterium avium and Mycobacterium intracellulare." J. Clin. Microbiol. 52 (2014): 4439-4440. PubMed: 25274991.
- Beggs, M. L., R. Stevanova and K. D. Eisenach. "Species Identification of *Mycobacterium avium* Complex Isolates by a Variety of Molecular Techniques." <u>J. Clin. Microbiol.</u> 38 (2000): 508-512. PubMed: 10655336.
- Devallois, A., et al. "Molecular Characterization of Mycobacterium avium Complex Isolates Giving Discordant Results in AccuProbe Tests by PCR- Restriction Enzyme Analysis, 16S rRNA Gene Sequencing, and DT1-DT6 PCR." J. Clin. Microbiol. 35 (1997): 2767-2772. PubMed: 9350730.
- Wayne, L. G., et al. "Serovar Determination and Molecular Taxonomic Correlation in Mycobacterium avium, Mycobacterium intracellulare, and Mycobacterium scrofulaceum: A Cooperative Study of the International Working Group on Mycobacterial Taxonomy." <u>Int. J. Syst.</u> Bacteriol. 43 (1993): 482-489. PubMed: 8347508.
- 7. Koh, W. J., et al. "Clinical Significance of the Differentiation between *Mycobacterium avium* and *Mycobacterium intracellulare* in *M. avium* Complex Lung Disease." Chest 142 (2012): 1482-1488. PubMed: 22628488.

ATCC<sup>®</sup> is a trademark of the American Type Culture Collection.

BEI Resources

www.beiresources.org

E-mail: <a href="mailto:contact@beiresources.org">contact@beiresources.org</a>

Tel: 800-359-7370

Tel: 800-359-7370 Fax: 703-365-2898